

	DOCUMENT TITLE: Installation Qualification Protocol, Weighing Scale	Revision: A
	DOCUMENT #: VV.SLQ017	Effective Date: 29Nov2022

PROTOCOL ATTACHMENT 1 – IQ: Verification of the installation

Verification	Acceptance criteria	Results observed	Conformity C / NC
Brand	Mettler Toledo	Mettler Toledo	C
Model number	AB54-S	AB54-S	C

Comments:	Equipment ID #: ST-012 Serial number: 1125311515
------------------	---

Completed By (Sign and Date): *Nalle* 5 Jan 2023

Reviewed By (Sign and Date): *Ethan Rao* 5 Jan 2023

	DOCUMENT TITLE: Installation Qualification Protocol, Weighing Scale	Revision: A
	DOCUMENT #: VV.SLQ017	Effective Date: 29Nov2022

PROTOCOL ATTACHMENT 2 – IQ: Verification of the documentation

Documentation	Acceptance criteria	Results observed	Conformity C / NC
Equipment user manual	Presence of the user manual in English	Manual is Present	C
Equipment calibration certificate	Presence of equipment calibration certificate	Calibration Certificate is Present	C
Calibration requirements	Schedule / frequency of calibration determined	1 Year frequency	C

Comments:	
------------------	--

Completed By (Sign and Date): *Nalle* 5 Jan 2023

Reviewed By (Sign and Date): *Ethan Rao* 5 Jan 2023

	DOCUMENT TITLE: Installation Qualification Protocol, Weighing Scale	Revision: A
	DOCUMENT #: VV.SLQ017	Effective Date: 29Nov2022

PROTOCOL ATTACHMENT 3 – IQ: Verification of utilities

Verification	Acceptance criteria	Results observed	Conformity C / NC
Power supply	Compliant if the weighing scale is powered by AC adapter cable	AC Adapter Cable is Present	C
Wiring	No damage to AC adapter cable	AC Adapter Cable is Undamaged	C
Operating temperature	10 – 30°C	22.9°C	C
Operating relative humidity	15 – 80% RH	50%	C
Operation	Equipment can be powered on and is operational	Equipment can be powered on and is operational	C

Comments:	
------------------	--

Completed By (Sign and Date): *Nalle* 5 Jan 2023

	DOCUMENT TITLE: Installation Qualification Protocol, Weighing Scale	Revision: A
	DOCUMENT #: VV.SLQ017	Effective Date: 29Nov2022

Reviewed By (Sign and Date): *Ethan Rao* 5 Jan 2023

	DOCUMENT TITLE: Installation Qualification Protocol, Weighing Scale	Revision: A
	DOCUMENT #: VV.SLQ017	Effective Date: 29Nov2022

PROTOCOL ATTACHMENT 4 – Calibration Log

In the table below, record the calibration information for all equipment used during the test. If the equipment is “calibrate before use” then include or refer to the procedure used for calibration, record any equipment required and its calibration info, and record the results of the calibration (data).

Instrument ID	Description	Use	Calibration Date	Due Date	Verified By/Date
ST-013	Temp/Humidity Monitor	Measure temp/humidity	10 Nov 2022	10 Nov 2024	Ethan Rao 5 Jan 2023

	DOCUMENT TITLE: Installation Qualification Protocol, Weighing Scale	Revision: A
	DOCUMENT #: VV.SLQ017	Effective Date: 29Nov2022

PROTOCOL ATTACHMENT 5 – Signature Log

In the table below, all personnel who have completed any section of this protocol must include the listed information.

Printed Name	Title	Initials	Signature	Company	Date
Ethan Rao	Director of R&D	ER	<i>Ethan Rao</i>	Silq Technologies	5 Jan 2023
Na He	Research Scientist	NH	<i>Na He</i>	Silq Technologies	5 Jan 2023