

Certificate of Conformance

Pathway Medtech LLC. certifies that the products or services indicated below have been manufactured in accordance with the requirements specified on referenced documents.

Part Number:	SLQ-211610SPT
Description:	SILQ 2-way Cleartract SPT 16 Fr x 10 ml Catheter, Sterile, 10-Pack
Revision:	B
Lot Number(s):	SLQ-05012025
Serial Number(s):	N/A
Expiration Date:	2028-05-31
Quantity Released:	1298 catheters/129.8 shippers


The product indicated can be transferred to a “release” inventory location and is cleared for distribution. The individual signing below authorizes the release of the product.


 Signature (Quality Assurance Representative)

02/15/2025
 Date

	Work Order Form	Doc No.: FM-0039
	This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.	Rev.: G
		Page 1 of 7



Production Planning – Document the products identification information with the planned build quantity and any deviation(s) that may apply. Attach a copy of the device's specification/drawing/BOM and any applicable deviation(s) to the work order form prior to transfer to manufacturing.

Part No.	Rev	Lot No./ Serial No.	Expiration Date	Quantity	Deviation No.
SLQ-211610SPT	B	SLQ-05012025	2028-05-31	2200	N/A
Comments					
SILQ 2-way Cleartract SPT 16 Fr x 10ml Catheter, Sterile, 10-Pack					
Name		Signature			Date
Sean Tamayo					05/14/2025

Lot No/Serial Number	SLQ-05012025
----------------------	--------------

	Work Order Form	Doc No.: FM-0039
	This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.	Rev.: G
		Page 2 of 7

Line Clearance - The individual performing the line clearance certifies that line clearance activities were conducted per WI-0014 prior to starting manufacturing activities. Quality management or designee verifies the line clearance was performed. *Note: The date line clearance occurred is considered the start date of manufacturing.*

Actions		Performed by	Date
1. Clear Area of previous WO Raw Materials.	☒ Completed.		05/14/2025
2. Verify kitted Raw Materials (RM) for: a. RM meet BOM revision level requirements. b. RM have not Expired. c. RM have no physical damage	a. ☒ Completed. b. ☒ Completed. c. ☒ Completed.		05/14/2025
3. Equipment is calibrated and properly functioning. <i>Note: Equipment identification and calibration due date will be documented in the work order.</i>	☒ Completed. ☐ N/A		05/14/2025
4. Proper documentation is present and available.	☒ Completed.		05/14/2025
5. Production Associates are trained and qualified to perform the process.	☒ Completed. ☐ N/A		05/14/2025
6. All work surfaces cleaned with IPA.	☒ Completed. ☐ N/A		05/14/2025
7. Labels printed per Label Generation and Reconciliation Form.	☐ Completed. ① ☒ N/A		05/14/2025
8. Verify validated processes have been set to critical points within specification.	☒ Completed. ☐ N/A	James Romero	05/14/2025
9. Verify Start-up activities are within specification	☒ Completed. ☐ N/A	James Romero	05/14/2025
Conclusion: Line clearance performed	☒ Completed		05/14/2025

Quality Management		
Name	Signature	Date
Giovanna Siqueros		07/15/2025

① Labels to be printed at a later date. ST 05/14/2025

Lot No/Serial Number	SLQ-05012025
----------------------	--------------

	Work Order Form	Doc No.: FM-0039
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.</i>	Rev.: G
		Page 3 of 7

Material/ Procedure Traceability (Check appropriate Box):

☐ – Identification information for the procedures and materials used correspond to the attached product's BOM and are captured on the attached Production Pick List (ref. P/N FM-0059).

☒ - Identification information for the procedures and materials used correspond to the attached product's BOM and are captured in the table below (N/A or delete table if not applicable).

BOM Item No.	Type	Rev	Material Designation	Lot No./Serial No.	Expiration Date	Rev	Material Designation	Lot No./Serial No.	Expiration Date
SLQ-8104	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	SLQ-05012025	2028-05-31		☐ Primary ☐ Substitute		
1	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	20220902	N/A		☐ Primary ☐ Substitute		
2	☐ Document ☒ Material	B	☒ Primary ☐ Substitute	2502-01-2	2027/02/12		☐ Primary ☐ Substitute		
3	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
4	☐ Document ☒ Material	A	☐ Primary ☒ Substitute	A199054	N/A		☐ Primary ☐ Substitute		
5	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	SLQ-05012025	2028-05-31		☐ Primary ☐ Substitute	N/A ST 06/19/2025	
6	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	000072 011223	N/A		☐ Primary ☐ Substitute		
7	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	SLQ-05012025	2028-05-31		☐ Primary ☐ Substitute		
8	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
9	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
10	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
11	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
12	☒ Document ☐ Material	B	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
13	☒ Document ☐ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
14	☒ Document ☐ Material	B	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		

Lot No/Serial Number	SLQ-05012025
----------------------	--------------



Work Order Form

Doc No.: FM-0039

Rev.: G

Page 4 of 7

This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.



BOM Item No.	Type	Rev	Material Designation	Lot No./Serial No.	Expiration Date	Rev	Material Designation	Lot No./Serial No.	Expiration Date
15	☒ Document ☐ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
16	☒ Document ☐ Material	B	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
SLQ-211610SPT	☐ Document ☒ Material	B	☒ Primary ☐ Substitute	SLQ-05012025	2028-05-31		☐ Primary ☐ Substitute		
1	☐ Document ☒ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
2	☒ Document ☐ Material	E	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
3	☒ Document ☐ Material	G	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute	N/A ST 06/19/2025	
4	☒ Document ☐ Material	B	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
5	☒ Document ☐ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
6	☒ Document ☐ Material	A	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
7	☒ Document ☐ Material	F	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
8	☒ Document ☐ Material	C	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
9	☒ Document ☐ Material	F	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		
10	☒ Document ☐ Material	B	☒ Primary ☐ Substitute	N/A	N/A		☐ Primary ☐ Substitute		

Lot No/Serial Number

SLQ-05012025

	Work Order Form		Doc No.: FM-0039
	This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.		Rev.: G
			Page 5 of 7

Training Traceability – All individuals indicated in the following table have been trained to the referenced document(s) which support the manufacturing, inspection and test activities of the product identified in the work order. Leadership representative has certified that the individual(s) have received proper training prior to conducting manufacturing activities.

Document No.	Rev	Print Name	Signature
MP-SLQ-0001	B	Saverdy Pen	
		Joseph Ragland	
		Brandy Mercado	
		Nicholas Spillman	
MP-SLQ-0002	B	Saverdy Pen	
		Joseph Ragland	
		Brandy Mercado	
		Nicholas Spillman	
Leadership Representative			
Name		Signature	Date
James Romero			05/17/2025

ST 05/22/2025

*Note, Brandy and Nicholas added to ~~the~~ the line
on 5/22/2025 ST 05/22/2025

Lot No/Serial Number	SLQ-05012025
----------------------	--------------

	Work Order Form		Doc No.: FM-0039
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.</i>		Rev.: G
			Page 6 of 7

Equipment Traceability – Identify any equipment that is used in the manufacturing/processing of the product in the following table.

Equipment Description	Asset No.	Cal Due Date (If applicable)	Equipment Description	Asset No.	Cal Due Date (if applicable)
Accu Seal Impulse Sealer	ES-0099-01	08/09/2025	Catheter Automation Machine	ES-0091-01	N/A

Production Reporting – Document scrapped information, date manufacturing completed, and final quantity manufactured.

BOM Item No.	Defect No. 1		Defect No. 2		Defect No. 3		Defect No. 4		Defect No. 5		Comments
	Code	QTY	Code	QTY	Code	QTY	Code	QTY	Code	QTY	
SLQ-6004	3	1140	10	48	N/A	N/A	N/A	N/A	N/A	N/A	Code 10 Breakdown: QTY:35 pulled for testing by Ethan Rao QTY:10 for Bioburden/LAL
N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Device Assembly	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

Defect Codes: 1=Damage, 2=Dimensional, 3=Visual, 4=Voids/Shorts, 5=Discoloration, 6=Contamination, 7=Operator Error, 8=Occlusion, 9=Leak, 10 =Miscellaneous (Indicate in Comment Section)

Manufacturing Completion Date	06/18/2025		Quantity MFG	1298 catheters/ 129.8 shippers
Name	Signature		Date	
Sean Tamayo			06/18/2025	

Lot No/Serial Number	SLQ-05012025
----------------------	--------------



Pathway Medtech, LLC.

Printed by Sean Tamayo

Printed on 05/08/2025

Local time zone (GMT-08:00) Pacific Time (America/Los Angeles)

Page Arena > Workspace Items > Bill of Materials > Indented

All information contained in this document is proprietary and confidential.

SLQ-211610SPT

In Production

SILQ 2-way Cleartract SPT 16 Fr x 10 ml Catheter, Sterile, 10-Pack


Revision »B - ECO-001669 Status: Effective, Status Date: 05/14/2024 03:24:42 PM, Shared

Contains **10** first-level Items, **50** line Items, **41** unique Items, **33** of which are shared.

# ITEM NUMBER	ITEM NAME	CATEGORY	PHASE	△	↑	□	↺	⬇	⬆	⬇	⬆	⬇	⬆
10 FM-0029 rev B	Form, Certificate of Conformance	Form	In Prod	↑	□	2	1 each						
3 FM-0039 rev G	Form, Work Order	Form	In Prod	↑	□	2	1 each						
5 MP-SLQ-0003 rev A	Manufacturing Process Flow, SILQ Foley Catheter	Manufacturing Procedure	In Prod	↑	□	2	1 EA						
9 QP-010 rev F	Final Product Inspection	SOP Quality System Procedure	In Prod	↑	□	6	1 each						
7 QP-032 rev F	EO Sterilization Process Control	SOP Quality System Procedure	In Prod	↑	□	8	1 each						
6 SLQ-0001 rev A	Sterilization Load Shipping Pallet Configuration	Subassembly	In Prod	↑	□	2	1 EA						
1 SLQ-8104 rev A	SILQ 2-way Cleartract SPT 16 Fr x 10 ml Catheter, Non-Sterile, 10-Pack	Subassembly	In Prod	↑	□	3	1 EA						
2 WI-0006 rev E	Work Order Process	Work Instruction	In Prod	↑	□	2	1 each						
4 WI-0014 rev B	Line Clearance Process	Work Instruction	In Prod	↑	□	2	1 each						
8 WI-0021 rev C	Work Instruction, Pyrogen (LAL) Evaluation Process	Work Instruction	In Prod	↑	□	2	1 each						

NOTES: UNLESS OTHERWISE SPECIFIED

1. GENERAL
 1. ANY CHANGES TO SPECIFICATIONS SHALL BE APPROVED BY SILQ PRIOR TO PRODUCTION.
 2. FINAL PRODUCT INSPECTION PER QP-010.
 3. PRODUCT EXPIRY DATE IS 36 MONTHS FROM THE DATE OF MANUFACTURE.
2. STERILIZATION
 1. PRODUCT IS STERILIZED VIA ETHYLENE OXIDE (ETO) STERILIZATION PER QP-032 CYCLE PWI-7374.
3. TESTING
 1. PYROGEN EVALUATION PER WI-0021 AT ≤ 0.5 EU/ML OR ≤ 20.0 EU/DEVICE



REVISIONS	
REV.	DESCRIPTION
A	INITIAL RELEASE
B	CHANGE PRODUCT EXPIRY DATE FROM 12 TO 36 MONTHS. REVISED FORMAT TO MATCH T-018

UNLESS OTHERWISE SPECIFIED:
 DIMENSIONS ARE IN MILLIMETERS
 TOLERANCES:
 ANGULAR: MACH ± 0.2 BEND ± 2.0
 NO PLACE DECIMAL ± 0.5
 ONE PLACE DECIMAL ± 0.3
 TWO PLACE DECIMAL ± 0.15
 THREE PLACE DECIMAL ± 0.05

PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PATHWAY MEDTECH. REPRODUCTION OF THIS DOCUMENT IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PATHWAY MEDTECH IS PROHIBITED.

pathway 

TITLE: 2-WAY CLEARTRACT CATHETER,
 16 FR X 10 ML, STERILE, 10-PACK

DRAWN BY:

B. SANDOVAL

SIZE DWG. NO.

B SLQ-211610SPT

REV

B

DRAWING NOT TO SCALE

T-018 Rev. A

SHEET 1 OF 1



Pathway Medtech, LLC.

Printed by Sean Tamayo
Printed on 05/09/2025
Local time zone (GMT-08:00) Pacific Time (America/Los Angeles)
Page Arena > Workspace Items > Bill of Materials > Indented
All information contained in this document is proprietary and confidential.

SLQ-8104 In Production

SILQ 2-way Cleartract SPT 16 Fr x 10 ml Catheter, Non-Sterile, 10-Pack

Revision »A - ECO-001669 Status: Effective, Status Date: 05/14/2024 03:24:42 PM, Shared

Contains 16 first-level Items, 31 line Items, 24 unique Items, 23 of which are shared.


# ITEM NUMBER	ITEM NAME	CATEGORY	PHASE			
1 SLQ-6004 rev A	Catheter, 16 Fr x 10 mL, Silicone, NY Foley, 2-Way	Fabricated Part	In Prod		3	10 EA
2 SLQ-3000 rev B	Solution, SILQ Polymer	Purchased Component	In Prod		4	80 ML
SUBSTITUTES						
	RANK ITEM NUMBER	ITEM NAME		NOTES		
	1 HDX-3000 rev A	Solution, Hydrophilix Polymer	80 ML			
3 SLQ-5000 rev A	Bag, 2 W x 18 L x 2 Mil Thk, Polyethylene	Packaging	In Prod		3	10 EA
4 SLQ-5001 rev A	Pouch, Chevron, 3.5 X 21, Amcor TLP-058	Packaging	In Prod		4	10 EA
SUBSTITUTES						
	RANK ITEM NUMBER	ITEM NAME		NOTES		
	1 HDX-5001 rev A	Pouch, Chevron, 3.5 X 21, Amcor TLP-058	10 EA			
▶ 5 SLQ-LB-4008 rev A	Label, Pouch, SILQ ClearTract SPT Catheter	Labeling	In Prod		2	10 EA

# ITEM NUMBER	ITEM NAME	CATEGORY	PHASE			
6 SLQ-5016 rev A	Shelf Carton, 3.25 x 2.5 x 19.75, .018 SBS, Pre-Printed	Packaging	In Prod		3	1 EA
SUBSTITUTES						
	RANK ITEM NUMBER	ITEM NAME		NOTES		
	1 SLQ-5015 rev A	Shelf Carton, 3.25 1 EA x 2.5 x 19.75, .018 SBS				
	2 SLQ-5014 rev A	Shelf Carton, 3.25 1 EA x 2.5 x 19.75, .024 SBS				
7 SLQ-LB-4009 rev A	Label, Box, SILQ ClearTract SPT Catheter	Labeling	In Prod		2	2 EA
8 SLQ-4007 rev A	Artwork, IFU, SILQ 2-Way ClearTract SPT Catheter	Labeling	In Prod		2	1 EA
9 SLQ-5012 rev A	Box, 20 in. L x 18 in. W x 6 in. D, 200 LB Test	Packaging	In Prod		2	0.1 EA
10 PWY-5000 rev A	Tape, Packaging, 2 Mil, 2 in x 110 yds, Clear	Packaging	In Prod		2	5.2 IN
11 SLQ-LB-4010 rev A	Label, Shipper Box, SILQ ClearTract SPT Catheter	Labeling	In Prod		2	0.2 EA
12 MP-SLQ-0002 rev B	Manufacturing Procedure, SILQ ClearTract Catheter Surface Treatment Using the Automation Machine	Manufacturing Procedure	In Prod		2	1 EA
13 FM-SLQ-0001 rev A	Form, Run Sheet for SILQ Automation Machine	Form	In Prod		2	1 EA
14 TM-HDX-0001 rev B	Test Method, UV Light Absorbance, Hydrophilix Catheters	Test Method	In Prod		2	1 each
15 TM-HDX-0002 rev A	Test Method, Rinse Test, Hydrophilix Catheter	Test Method	In Prod		2	1 each
16 MP-SLQ-0001 rev B	Manufacturing Procedure, Labeling and Packaging, SILQ ClearTract Catheters	Manufacturing Procedure	In Prod		2	1 EA

NOTES: UNLESS OTHERWISE SPECIFIED

1. GENERAL
 1. ANY CHANGES TO SPECIFICATIONS SHALL BE APPROVED BY PATHWAY PRIOR TO PRODUCTION.
 2. DRAWING IS FOR REFERENCE ONLY. REFER TO THE RELEASED BILL OF MATERIALS FOR THE COMPLETE PARTS LIST.
 3. UP UNTIL SEALING THE POUCH, PRODUCT SHALL BE PACKAGED IN A CLEANROOM ENVIRONMENT. POST SEALING THE POUCH, ASSEMBLY WILL CONTINUE OUTSIDE OF THE CLEANROOM ENVIRONMENT.
2. ASSEMBLY PROCEDURE
 1. APPLY POLYMER SOLUTION TO THE CATHETER IN THE HATCHED AREA SHOWN. ALLOW THE CATHETER TO DRY OVERNIGHT.
 2. INSERT THE CATHETER INTO A POLYETHYLENE SLEEVE SUCH THAT THE TIP OF THE CATHETER IS ALIGNED TOWARDS THE CLOSED END.
 3. PLACE THE POUCH LABEL ON THE CLEAR SIDE OF THE POUCH. WHEN VIEWING THE CLEAR SIDE OF THE POUCH FROM ABOVE, WITH THE CHEVRON TO THE RIGHT, THE LABEL TEXT SHALL READ LEFT TO RIGHT. WITH BEST EFFORT, APPLY THE LABEL SUCH THAT IT IS CENTERED.
 4. INSERT THE SLEEVE / CATHETER INTO THE POUCH SUCH THAT THE TIP OF THE CATHETER IS ALIGNED TOWARDS THE CHEVRON END OF THE POUCH.
 5. SEAL THE POUCH.
 6. PLACE 10 SLEEVED, POUCHED, SEALED, AND LABELED CATHETERS INTO ONE SHELF CARTON. APPLY TWO SHIPPER BOX LABELS TO THE LONG, WIDER SIDE OF THE CARTON. VERIFY THAT THE TEXT READS LEFT TO RIGHT WHEN THE OPENING OF THE CARTON IS TO THE LEFT. WITH BEST EFFORT, CENTER THE LABEL ON ITS RESPECTIVE FACE OF THE CARTON.
 7. PLACE 10 SHELF CARTONS INTO ONE SHIPPER BOX. SEAL WITH PACKAGING TAPE AS NEEDED. PLACE A SHIPPER BOX LABEL IN THE UPPER LEFT CORNER ON TWO ADJACENT FACES, READING LEFT TO RIGHT.
3. DIMENSIONAL
 1. () INDICATES LABEL REFERENCE DIMENSIONS.
 2.  INDICATES CRITICAL DIMENSIONS.
4. CLEANING AND HANDLING
 1. ALL PRODUCT, UP UNTIL INSERTING THE POUCHES INTO SHELF CARTONS, SHALL BE HANDLED BY GLOVED HANDS APPROPRIATE FOR A CLEANROOM ENVIRONMENT.
5. FUNCTIONAL INSPECTION (AQL = 6.5 UNLESS OTHERWISE STATED)
 1. THE CATHETER'S ABSORBANCE VALUE AT 585 NANOMETERS SHALL BE GREATER THAN OR EQUAL TO 0.300 AU. REFER TO TM-HDX-0001.
 2. THE CATHETER'S ABSORBANCE VALUE AT 440 NANOMETERS SHALL BE LESS THAN OR EQUAL TO 0.020 AU. REFER TO TM-HDX-0002.
 3. THE PEEL STRENGTH FOR A 1 INCH WIDE SEALED SECTION OF THE STERILE BARRIER PACKAGING SHALL BE ≥ 1.00 LB WHEN TESTED PER TM-0003. IT IS NOT REQUIRED TO PLACE CATHETERS IN THESE TEST POUCHES, OR LABEL THEM. INSPECT N=3 SAMPLES, 1 AT THE BEGINNING, MIDDLE, AND END OF THE MANUFACTURING RUN.
6. VISUAL INSPECTION (AQL = 6.5 UNLESS OTHERWISE STATED)
 1. THE CATHETER SHALL BE FREE OF CURED POLYMER AND DEBRIS.
 2. THE CATHETER SHALL BE PROPERLY ORIENTED IN THE SLEEVE AND POUCH WITH THE LABEL APPLIED TO THE APPROPRIATE FACE AND IN THE CORRECT ORIENTATION.
 3. INSIDE THE STERILE BARRIER, DEVICE IS TO BE FREE OF PARTICULATES WHEN VIEWED AT 18 INCHES AWAY UNDER AMBIENT ROOM LIGHTING AND WITHOUT THE USE OF MAGNIFICATION.
 4. THE STERILE BARRIER PACKAGING SHALL NOT INCLUDE ANY OF THE FOLLOWING DEFECTS:
 1. UNSEALED AREAS, NONHOMOGENEOUS OR UNSEALED AREAS, OVERSEALED AREAS, ANROW SEALS, CHANNELS, WRINKLES / FOLDOVERS / CRACKS, OR TEARS / PINHOLES.
 5. THE SHELF CARTON SHALL HAVE THE CORRECT NUMBER OF CATHETERS INSIDE AND FACING THE CORRECT ORIENTATION. THE LABELS SHALL BE AFFIXED IN THE CORRECT LOCATION AND ORIENTATION.
 6. THE SHIPPER BOX SHALL HAVE THE CORRECT NUMBER OF SHELF CARTONS. THE LABELS SHALL BE APPLIED IN THE CORRECT LOCATION AND ORIENTATION.

REVISIONS			
ZONE	REV.	DESCRIPTION	DATE
	A	INITIAL RELEASE	



UNLESS OTHERWISE SPECIFIED:
 DIMENSIONS ARE IN MILLIMETERS
 TOLERANCES
 ANGULAR: MACH ± 0.2 BEND ± 0.0
 NO PLACE DECIMAL ± 0.5
 ONE PLACE DECIMAL ± 0.3
 TWO PLACE DECIMAL ± 0.15
 THREE PLACE DECIMAL ± 0.05

PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PATHWAY NPI. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PATHWAY NPI IS PROHIBITED.


pathway
 New Product Introduction

TITLE:
 CATHETER, SILICONE, NY FOLEY, NON-STERILE, 10-PACK, MASTER DRAWING

SIZE DWG. NO. REV
B SLQ-81XX A

DO NOT SCALE DRAWING SCALE: 1:5 SHEET 1 OF 3





ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	SLQ-5001	POUCH, CHEVRON, 3.5 X 21, AMCOR TLP-058	10
2	SLQ-5000	BAG, 2 W X 18 L X 2 MIL THK, POLYETHYLENE	10
3	SLQ-60XX	SEE TABLE FOR PART SPECIFICATION	10
4	SLQ-5016	SHELF CARTON, 3.25 X 2.5 X 19.75, .018 SBS, PRE-PRINTED	1
5	SLQ-LB-4008	POUCH LABEL, CLEARTRACT SP CATHETER	10
6	SLQ-LB-4009	SHELF CARTON LABEL, CLEARTRACT SP CATHETER	1
7	SLQ-4007	IFU, SILQ 2-WAY CLEARTRACT SP CATHETER	1



PART NUMBER	DESCRIPTION	SLQ-81XX
SLQ-6000	CATHETER, 14 FR X 10 ML, SILICONE, NY FOLEY, 2-WAY	SLQ-8100
SLQ-6001	CATHETER, 16 FR X 5 ML, SILICONE, NY FOLEY, 2-WAY	SLQ-8101
SLQ-6002	CATHETER, 18 FR X 10 ML, SILICONE, NY FOLEY, 2-WAY	SLQ-8102
SLQ-6004	CATHETER, 16 FR X 10 ML, SILICONE, NY FOLEY, 2-WAY	SLQ-8104

FDA
REGISTERED

ISO
CERTIFIED 13485

SIZE	DWG. NO.	REV
B	SLQ-81XX	A
SCALE: 1:1		SHEET 2 OF 3

ITEM NO.	PARTNO	DESCRIPTION	QTY.
8	SLQ-5012	BOX, 20 IN. L X 18 IN. W X 6 IN. D, 200 LB TEST	1
9	PWY-5000	TAPE, PACKAGING, 2 MIL, 2 IN X 110 YDS, CLEAR	AR
10	SLQ-LB-4010	SHIPPER BOX LABEL, CLEARTRACT SP CATHETER	2



FDA
REGISTERED

ISO
CERTIFIED 13485

SIZE	DWG. NO.	REV
B	SLQ-81XX	A
SCALE: 1:1		SHEET 3 OF 3

	Work Order Form	Doc No.: FM-0039
	This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.	Rev.: G
		Page 7 of 7

Documentation Checklist – If applicable, the following documentation shall be included upon completion of the manufacturing build. Document packet shall be organized in the order indicated below. Ensure that good documentation practices are applied per QP-003 "Quality Records". Check box(es) ☐ that apply to each requirement. The individual who signs has confirmed that all requirements have been met.

1. Work Order (WO) Record	☒ Present and Completed	
2. Product's Specification/Drawing & BOM (Arena PLM)	☒ Present	
3. Deviations	☐ Present	☒ N/A
4. Product Labeling & Generation/ Reconciliation Record	☒ Present, Completed and Results Meet Specification	☐ N/A
5. Non-Conforming Reports (NCR)	☐ Present and Completed	☒ N/A
6. Inspection Records	☒ Present, Completed and Results Meet Specification	
7. Critical Process Parameter Record	☒ Present, Completed and Results Meet Specification	☐ N/A
Name	Signature	Date
Giovanna Siqueiros		07/15/2025

Release Authorization for Further Processing (e.g., terminal sterilization) – The individual signing below is certifying that they independently reviewed the Device History Record (DHR) and confirmed that all required documentation is included, records are complete, good documentation practices were followed, and specifications have been met (if not NCR is included)

Name	Signature	Date
Giovanna Siqueiros		07/15/2025

Lot No/Serial Number	SLQ-05012025
----------------------	--------------



Form, Run Sheet for SILQ Automation Machine

This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from. If any portion of this document is in hardcopy form, it must be considered and treated as an uncontrolled "For Reference Only", unless accompanied with a controlled stamp. Users of this document must verify that they are using the most current revision of the controlled document before performing work.

Doc No.: FM-SQL-0001

Rev.: A

Page 1 of 1


Part Number	Revision	Description	Machine Asset No.
SLQ-8104	A	SILQ 2-Way ClearTract SPT 16Fr x 10ml Catheter, Non-Sterile, 10-Pack	ES-0091-01

Equipment Setup / Startup		Adjustments during Production			
Setting [Units]		Value	Value	Value	Value
Engineering / Lead Supervisor	Setup / Adjustment Date	05/14/2025	N/A	ST 06/18/2025	
	Setup / Adjustment Time	10:45			
	Fill Time [seconds]	7.2			
	Confirm Cylinder Actuation All 4 pairs [P/F]	P			
	UV Cure Time [seconds]	120			
	Drain Time [seconds]	30			
	Rinse Time [seconds]	3.8			
	Final Drain [seconds]	30			
	Number of Drain / Rinse Cycles [Count]	3			
	Confirm Cylinder Actuation All 4 pairs [P/F]	P			
	Size Setting	16 FR			
	Performed By Name	James Romero			
Performed By Signature	James Romero				
QA	Verified By Name	Giovanna Siqueiros			
	Verified By Signature	C. Siqueiros			

Date Inspected		Part No.		Rev.	Description				
05/14/2025-06/18/2025		SLQ-8104		A	SILQ 2-Way Cleartract SPT 16 Fr x 10 mL Catheter, Non-Sterile, 10-Pack				
Lot No.		Lot QTY		Approved Supplier(s)					
SLQ-05012025		1298		Pathway					
In-Process Inspection	Catheter has been coated with SILQ's polymer solution (SLQ-3000)	Catheter has been rinsed of cured polymer solution and is free of debris	Coated catheter absorbance @585nm ≥0.300 AU Ref TM-HDX-0001	Rinsed catheter absorbance @440nm ≤0.020 AU Ref TM-HDX-0002	Pouch Label Inspection: 1) Label content matches DMR 2) P/N, Lot No., QTY 1, and Exp date on labels match Work Order requirements. 3) Labels are legible, not missing text, 4) Labels with no stains or discoloration. 5) Non HRI Text exp date format is "YYYY-MM-DD".	Carton Label Inspection: 1) Label content matches DMR 2) P/N, Lot No., QTY 10, and Exp date on labels match Work Order requirements. 3) Labels are legible, not missing text, 4) Labels with no stains or discoloration. 5) Non HRI Text exp date format is "YYYY-MM-DD".	Pouch Inspection: 1) Catheter is properly oriented in the sleeve and pouch. 2) Pouch seal is continuous with no wrinkles, gaps, folds, bubbles or burn/overheating clear spots. Seal width at least 1/8 inch. 3) Pouch is with no tears, punctures or pin holes. 4) Pouch observed with no particulate/loose debris inside pouch. 5) Labels affixed & positioned in correct location, is legible.	Inner Carton Inspection: 1) Correct quantity in 10-Pack Box. 2) Product configuration matches specification. 3) Labels affixed & positioned in correct location. 4) No stains or discoloration that affects appearance. 5) No damage observed to the device, its sterile barrier, or inner carton.	Shipper Box Inspection: 1) Correct quantity in Box. 2) Product configuration matches specification. 3) Labels affixed & positioned in correct location. 4) No stains or discoloration that affects appearance. 5) No damage observed to the inner carton, or shipper box. IP-0002
	Revision	A	A	B	A	A	A	A	B
	Lower Spec Limit (LSL)	N/A	N/A	0.300	N/A	N/A	N/A	N/A	N/A
	Upper Spec Limit (USL)	N/A	N/A	N/A	0.020	N/A	N/A	N/A	N/A
	UOM	EA	EA	AU	AU	EA	EA	EA	EA
	Measuring Equipment Type	Visual	Visual	Spectrophotometer	Spectrophotometer	Visual	Visual	Visual	Visual
	Asset No.	N/A	N/A	See Comments	See Comments	N/A	N/A	N/A	N/A
Cal Due Date	N/A	N/A	See Comments	See Comments	N/A	N/A	N/A	N/A	
AQL	6.5	6.5	6.5	6.5	6.5	6.5	6.5	6.5	
Sample Size	13	13	N/A	N/A	13	13	13	13	
Sample No.		Variable/ Attribute Data							
1	P	P	N/A	N/A	P	P	P	P	P
2	P	P	N/A	N/A	P	P	P	P	P
3	P	P	N/A	N/A	P	P	P	P	P
4	P	P	N/A	N/A	P	P	P	P	P
5	P	P	N/A	N/A	P	P	P	P	P
6	P	P	N/A	N/A	P	P	P	P	P
7	P	P	N/A	N/A	P	P	P	P	P

8	P	P	N/A	N/A	P	P	P	P	P
9	P	P	N/A	N/A	P	P	P	P	P
10	P	P	N/A	N/A	P	P	P	P	P
11	P	P	N/A	N/A	P	P	P	P	P
12	P	P	N/A	N/A	P	P	P	P	P
13	P	P	N/A	N/A	P	P	P	P	P
Result Summary									
Pass QTY	13	13	N/A	N/A	13	13	13	13	13
Fail QTY	0	0	N/A	N/A	0	0	0	0	0
Result (P/F)	P	P	N/A	N/A	P	P	P	P	P
Statistical Summary									
Min	N/A	N/A	0.000	0.000	N/A	N/A	N/A	N/A	N/A
Max	N/A	N/A	0.000	0.000	N/A	N/A	N/A	N/A	N/A
Mean	N/A	N/A	#DIV/0!	#DIV/0!	N/A	N/A	N/A	N/A	N/A
StDev	N/A	N/A	#DIV/0!	#DIV/0!	N/A	N/A	N/A	N/A	N/A
Comments									
Rinse test and UV Absorbance Test were performed at SILQ. The rest of the visual inspections were performed at Pathway.									
Disposition:	ACCEPT					NCMR No:	N/A		
Inspected By:	Brandy Mercado					Date:	06/18/2025		
Verified By:	[Signature]					Date:	06/18/2025		

Date Inspected		Part No.		Rev.	Description				
6/13/2025		SLQ-8104		A	SILQ 2-Way Cleartract SPT 16 Fr x 10 mL Catheter, Non-Sterile, 10-Pack				
Lot No.		Lot QTY		Approved Supplier(s)					
SLQ-05012025		1311 1298		ST 06/18/2025 Pathway					
In-Process Inspection	Catheter has been coated with SILQ's polymer solution (SLQ-3000)	Catheter has been rinsed of cured polymer solution and is free of debris	Coated catheter absorbance @585nm ≥0.300 AU Ref TM-HDX-0001	Rinsed catheter absorbance @440nm ≤0.020 AU Ref TM-HDX-0002	Pouch Label Inspection: 1) Label content matches DMR 2) P/N, Lot No., QTY 1, and Exp date on labels match Work Order requirements. 3) Labels are legible, not missing text, 4) Labels with no stains or discoloration. 5) Non HRI Text exp date format is " YYYY-MM-DD ".	Carton Label Inspection: 1) Label content matches DMR 2) P/N, Lot No., QTY 10, and Exp date on labels match Work Order requirements. 3) Labels are legible, not missing text, 4) Labels with no stains or discoloration. 5) Non HRI Text exp date format is " YYYY-MM-DD ".	Pouch Inspection: 1) Catheter is properly oriented in the sleeve and pouch. 2) Pouch seal is continuous with no wrinkles, gaps, folds, bubbles or burn/overheating clear spots. Seal width at least 1/8 inch. 3) Pouch is with no tears, punctures or pin holes. 4) Pouch observed with no particulate/loose debris inside pouch. 5) Labels affixed & positioned in correct location, is legible.	Inner Carton Inspection: 1) Correct quantity in 10-Pack Box. 2) Product configuration matches specification. 3) Labels affixed & positioned in correct location. 4) No stains or discoloration that affects appearance. 5) No damage observed to the device, its sterile barrier, or inner carton.	Shipper Box Inspection: 1) Correct quantity in Box. 2) Product configuration matches specification. 3) Labels affixed & positioned in correct location. 4) No stains or discoloration that affects appearance. 5) No damage observed to the inner carton, or shipper box. IP-0002
	Revision	A	A	B	A	A	A	A	B
	Lower Spec Limit (LSL)	N/A	N/A	0.300	N/A	N/A	N/A	N/A	N/A
	Upper Spec Limit (USL)	N/A	N/A	N/A	0.020	N/A	N/A	N/A	N/A
	UOM	EA	EA	AU	AU	EA	EA	EA	EA
	Measuring Equipment Type	Visual	Visual	Spectrophotometer	Spectrophotometer	Visual	Visual	Visual	Visual
	Asset No.	N/A	N/A	See Comments	See Comments	N/A	N/A	N/A	N/A
Cal Due Date	N/A	N/A	See Comments	See Comments	N/A	N/A	N/A	N/A	
AQL	6.5	6.5	6.5	6.5	6.5	6.5	6.5	6.5	
Sample Size	13	13	13	13	13	13	13	13	
Sample No.		Variable/ Attribute Data							
1	N/A	N/A	0.434	0.006	N/A	N/A	N/A	N/A	N/A
2	N/A	N/A	0.461	0.003	N/A	N/A	N/A	N/A	N/A
3	N/A	N/A	0.457	0.007	N/A	N/A	N/A	N/A	N/A
4	N/A	N/A	0.498	0.001	N/A	N/A	N/A	N/A	N/A
5	N/A	N/A	0.435	0.005	N/A	N/A	N/A	N/A	N/A
6	N/A	N/A	0.465	0.000	N/A	N/A	N/A	N/A	N/A
7	N/A	N/A	0.487	0.008	N/A	N/A	N/A	N/A	N/A
8	N/A	N/A	0.438	0.006	N/A	N/A	N/A	N/A	N/A
9	N/A	N/A	0.461	0.001	N/A	N/A	N/A	N/A	N/A

10	N/A	N/A	0.444	0.009	N/A	N/A	N/A	N/A	N/A
11	N/A	N/A	0.449	0.010	N/A	N/A	N/A	N/A	N/A
12	N/A	N/A	0.482	0.011	N/A	N/A	N/A	N/A	N/A
13	N/A	N/A	0.468	0.006	N/A	N/A	N/A	N/A	N/A
Result Summary									
Pass QTY	N/A	N/A	13	13	N/A	N/A	N/A	N/A	N/A
Fail QTY	N/A	N/A	0	0	N/A	N/A	N/A	N/A	N/A
Result (P/F)	N/A	N/A	P	Pathway	P	N/A	N/A	N/A	N/A
Statistical Summary									
Min	N/A	N/A	0.434	0.000	N/A	N/A	N/A	N/A	N/A
Max	N/A	N/A	0.498	0.011	N/A	N/A	N/A	N/A	N/A
Mean	N/A	N/A	0.460	0.006	N/A	N/A	N/A	N/A	N/A
StDev	N/A	N/A	0.020	0.004	N/A	N/A	N/A	N/A	N/A
Comments									
Rinse test and UV Absorbance Test were performed at SILQ. The remaining inspections were performed at Pathway.									
Disposition:	ACCEPT					NCMR No:	N/A		
Inspected By:	Ethan Rao					Date:	13-Jun-25		
Verified By:						Date:	06/18/2025		

Date Inspected	Part No.	Rev.	Description						
06/04/2025-06/05/2025	SLQ-8104	A	SILQ 2-Way Cleartract SPT 16 Fr x 10 mL Catheter, Non-Sterile, 10-Pack						
Lot No.	Lot QTY	ES-0059-01 Emplex Model MPS6300-M Continuous Band Sealer - Validated Range: Temperature = 142 °C - 144 °C, Speed = 175 in/min., Pressure = 89 psi ES-0099-01 ACCU-SEAL Model 5300-15-B, Impulse Sealer - Validated Range: Temperature = 270 °F - 280 °F, Dwell = 4 sec., Pressure = 55 psi							
SLQ-05012025	1298								
Specification	TM-0003								
Revision	C								
PRODUCT PREP FOR TEST									
* Seal one sample pouch P/N SLQ-5001 at the start, middle and end of each shift.									
*Cut two 1.00 inch wide x 2.00 inch long strips from the manufacturer sealed end of the pouch.									
* Execute peel strength test per TM-0003 and record results.									
Specification	≥1.0	≥1.0	142-144 °C	175 in/min	89 psi	270 - 280 °F	4.0	55.000	
Lower Spec Limit (LSL)	1.0	1.0	142.0	175.0	89	270.0	4.0	55	
Upper Spec Limit (USL)	N/A	N/A	144.0	175.0	89	280.0	4.0	55	
UOM	lb	lb	°C	in/min	psi	°F	seconds	psi	
Measuring Equipment Type	Force Gauge	Force Gauge	Visual BandSealer Gauge ONLY ES-0059-01	Visual BandSealer Gauge ONLY ES-0059-01	Visual BandSealer Gauge ONLY ES-0059-01	Visual ImpulseSealer Gauge ONLY ES-0099-01	Visual ImpulseSealer Gauge ONLY ES-0099-01	Visual ImpulseSealer Gauge ONLY ES-0099-01	
Asset No.	ES-0028-01	ES-0028-01	N/A	N/A	N/A	ES-0099-01	ES-0099-01	ES-0099-01	
Cal Due Date	12/28/2024	12/28/2024	N/A	N/A	N/A	8/9/2025	8/9/2025	8/9/2025	
Sample ID	Top Strip	Bottom Strip	N/A	N/A	N/A	1	1	1	
Sample No.	Variable/ Attribute Data								
1: Start (Top Strip)	4.908	N/A	N/A	N/A	N/A	270.0	4.0	55.0	
2: Start (Bottom Strip)	N/A	5.210	N/A	N/A	N/A	N/A	N/A	N/A	
3: Middle (Top Strip)	6.022	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
4: Middle (Bottom Strip)	N/A	5.550	N/A	N/A	N/A	N/A	N/A	N/A	
5: End (Top Strip)	4.838	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
6: End (Bottom Strip)	N/A	6.976	N/A	N/A	N/A	N/A	N/A	N/A	
1: Start (Top Strip)	5.466	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
2: Start (Bottom Strip)	N/A	6.324	N/A	N/A	N/A	N/A	N/A	N/A	
3: Middle (Top Strip)	5.002	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
4: Middle (Bottom Strip)	N/A	4.984	N/A	N/A	N/A	N/A	N/A	N/A	
5: End (Top Strip)	6.534	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
6: End (Bottom Strip)	N/A	4.950	N/A	N/A	N/A	N/A	N/A	N/A	
Result Summary									
Pass QTY	6	6	N/A	N/A	N/A	1	1	1	
Fail QTY	0	0	N/A	N/A	N/A	0	0	0	
Result (P/F)	P	P	N/A	N/A	N/A	P	P	P	
Statistical Summary									
Min	4.838	4.950	#REF!	#REF!	#REF!	#REF!	#REF!	#REF!	
Max	6.534	6.976	#REF!	#REF!	0.00	#REF!	#REF!	#REF!	
Mean	5.462	5.666	#REF!	#REF!	0.00	#REF!	#REF!	#REF!	
StDev	0.689	0.818	#REF!	#REF!	#REF!	#REF!	#REF!	#REF!	
Comments									
N/A									
Disposition:	ACCEPT					NCMR No:	N/A		
Inspected By:						Date:	06/05/2025		

	Label Generation & Reconciliation				Doc No.: FM-0033	
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.</i>				Rev.: D	
					Page 1 of 2	

Label Information						
Part No.	SLQ-LB-4009	Rev.	A	Lot No.	SLQ-05012025 ☐ N/A	Label QTY
						10 ☐ N/A
Description	Label, Box, SILQ ClearTract SPT Catheter	UDI No.	10850041415151 ☐ N/A		Expiration Date:	2028-05-31 ☐ N/A

Label Quantity Calculation	
Required Label Quantity (A)	135
2% Overage (A×0.02=B)	3
QC Retain Labels	2
Quantity of Labels Required (A+B+2=C)	140
Quantity of Labels Printed (C=D)	140
Production's Initial	ST
Date	06/03/2025

Label Verification	
What is the revision of the label per the BOM?	A
What is the revision on the printed label?	A
Is the printed information correct per specifications?	☒ Yes ☐ No
First and last label printed is affixed to back of form	☒ Yes ☐ No
Quantity of labels released to production (E)	138
QC's Initials	AS
Date	06/4/2025

Label Reconciliation	
Quantity of labels used by production (F)	135
Quantity scrapped by production (G)	3
Quantity of unused labels returned to QC and scrapped (H)	138
(E-F-G-H=0)	0
Are all labels accounted for?	☒ Yes ☐ No

If the counts cannot be reconciled, perform the following activities:

1. Verify that all products are labeled with the information recorded in the "Label Generation" section. Document the product quantity verified in designated field below.
2. Review the appropriate workstations and surrounding areas for the missing labels.

Investigational Results
☒ N/A

	Label Generation & Reconciliation This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.	Doc No.: FM-0033 Rev.: D Page 2 of 2
	QC's Initials <i>AS</i>	Date <i>07/15/2025</i>

QUANTITY	SIZE
10	16 Fr, 10 mL

2-WAY 100% SILICONE
CLEARTRACT SPT CATHETER



Rx Only ~~XX~~  



(01) 10850041415151
(10) SLQ-05012025
(17) 280531

Symbols are defined in Instruction for Use (IFU)
Store away from direct light exposure at room temperature
Not made with natural rubber latex



SILQ Technologies Corp.
323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

STERILE EO



MR

LOT

SLQ-05012025

REF

211610SPT



2028-05-31

SLQ-LB-4009 Rev. A

First

QUANTITY	SIZE
10	16 Fr, 10 mL

2-WAY 100% SILICONE
CLEARTRACT SPT CATHETER



Rx Only ~~XX~~  



(01) 10850041415151
(10) SLQ-05012025
(17) 280531

Symbols are defined in Instruction for Use (IFU)
Store away from direct light exposure at room temperature
Not made with natural rubber latex



SILQ Technologies Corp.
323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

STERILE EO



MR

LOT

SLQ-05012025

REF

211610SPT



2028-05-31

SLQ-LB-4009 Rev. A


Last

NOTES: UNLESS OTHERWISE SPECIFIED

1. GENERAL
 1. PRIOR TO PRINTING, VERIFY THE INFORMATION ON THE LABEL IS ACCURATE WHEN COMPARED AGAINST THE WORK ORDER.
 2. PRIOR TO PRINTING, VERIFY THE LABEL HAS GONE THROUGH THE INSTALLATION QUALIFICATION (IQ) PROCESS FOR LABELS. REFER TO FM-0044.
 3. LABEL CONTENT SHALL BE DEVELOPED PER WI-0009.
2. MATERIAL
 1. LABEL: STANDARD WHITE THERMAL TRANSFER
 2. RIBBON: WAX / RESIN
3. DIMENSIONS
 1. () INDICATES LABEL REFERENCE DIMENSIONS.
4. SAMPLING PLAN
 1. IN ADDITION TO THE LABEL GENERATION AND RECONCILIATION FORM, FM-0033, INSPECT N=3 ADDITIONAL SAMPLES AT THE BEGINNING OF THE RUN.
5. GS1 DATA MATRIX VALUES
 1. (01): GLOBAL TRADE ITEM NUMBER (GTIN)
 2. (10): BATCH OR LOT NUMBER
 3. (17): EXPIRATION DATE
6. VISUAL INSPECTION - LABEL TEXT
 1. POST-PRINTING, ALL TEXT ON THE LABEL SHALL BE LEGIBLE.
 2. ALL VARIABLE LABEL CONTENT (LOT NUMBER, EXPIRATION DATE, ETC.) SHALL MATCH THE WORK ORDER.
 3. THE GS1 DATA MATRIX SHALL BE FREE OF MISSING PIXELS AND SHALL NOT HAVE SMUDGES OR STREAKS THROUGH IT. WHEN SCANNED WITH A BARCODE SCANNER, THE TEXT ADJACENT TO THE GS1 DATA MATRIX SHALL MATCH THE TEXT DISPLAYED BY THE SCANNER.
 4. (01): VERIFY THE GTIN IS DETAILED IN A HUMAN READABLE FORMAT ADJACENT TO THE GS1 DATA MATRIX. THE GTIN SHALL BE A 14 DIGIT NUMBER.
 5. (10): VERIFY THE BATCH OR LOT NUMBER FORMAT IS IN COMPLIANCE WITH QP-019.
 6. (17): VERIFY THE EXPIRATION DATE IS IN THE FOLLOWING FORMAT: YYMMDD.
 7. VERIFY THE CATHETER STYLE AND CATHETER SIZE ARE REPRESENTATIVE OF A CONFIGURATION FOUND IN THE SKU TABLE.
 8. VERIFY THE GTIN IS FOUND IN THE SKU TABLE AND CORRESPONDS WITH THE CORRECT CATHETER SIZE AND STYLE.
 9. VERIFY THE REFERENCE NUMBER FOR THE CATHETER CORRESPONDS TO THAT FOUND IN THE SKU TABLE.
 10. VERIFY "QUANTITY" IS 10.
7. VISUAL INSPECTION - LABEL SYMBOLS
 1. POST-PRINTING, ALL SYMBOLS ON THE LABEL SHALL BE FREE OF SMUDGES / MISSING ARTWORK.
 2. THE EXPIRATION DATE NEXT TO ITS CORRESPONDING SYMBOL SHALL BE IN THE FORMAT: YYYY-MM-DD.

REVISIONS			
ZONE	REV.	DESCRIPTION	DATE
	A	INITIAL RELEASE	

ITEM NO.	PARTNO	DESCRIPTION	QTY.
1	SLQ-AW-4009	LABEL, BOX, SILQ CLEARTRACT SPT CATHETER	1
2	PWY-5007	RIBBON, THERMAL TRANSFER, WAX/RESIN, 4.3W	1
3	SLQ-5004	LABEL, 2.3125 IN. X 7.50 IN.	1



UNLESS OTHERWISE SPECIFIED:
 DIMENSIONS ARE IN MILLIMETERS
 TOLERANCES:
 ANGULAR: MACH ±0.2 BEND ±2.0
 NO PLACE DECIMAL ±0.5
 ONE PLACE DECIMAL ±0.3
 TWO PLACE DECIMAL ±0.15
 THREE PLACE DECIMAL ±0.05

PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PATHWAY NPL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PATHWAY NPL IS PROHIBITED.



REGISTERED CERTIFIED 13485


DO NOT SCALE DRAWING

pathway
 New Product Introduction

TITLE:
 LABEL, BOX, SILQ CLEARTRACT SPT CATHETER

SIZE DWG. NO. REV
B SLQ-LB-4009 A

SCALE: 1:1 SHEET 1 OF 2



SIZE	DWG. NO.	REV
B	SLQ-LB-4009	A
SCALE: 1:2		SHEET 2 OF 2

	Label Generation & Reconciliation		Doc No.: FM-0033
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.</i>		Rev.: D
			Page 1 of 2

Label Information					
Part No.	SLQ-LB-4008	Rev.	A	Lot No.	SLQ-05012025 ☐ N/A
Description	Label, Pouch, SILQ ClearTract SPT Catheter	UDI No.	00850041415154 ☐ N/A		Expiration Date:
				2028-05-31 ☐ N/A	Label QTY
				1 ☐ N/A	

Label Quantity Calculation	
Required Label Quantity (A)	1344
2% Overage (A×0.02=B)	27
QC Retain Labels	2
Quantity of Labels Required (A+B+2=C)	1373
Quantity of Labels Printed (C=D)	1373
Production's Initial	ST
Date	06/02/2025

Label Verification	
What is the revision of the label per the BOM?	A
What is the revision on the printed label?	A
Is the printed information correct per specifications?	☒ Yes ☐ No
First and last label printed is affixed to back of form	☒ Yes ☐ No
Quantity of labels released to production (E)	1371
QC's Initials	BS
Date	06/3/2025

Label Reconciliation	
Quantity of labels used by production (F)	1344
Quantity scrapped by production (G)	27
Quantity of unused labels returned to QC and scrapped (H)	0
(E-F-G-H=0)	0
Are all labels accounted for?	☒ Yes ☐ No

If the counts cannot be reconciled, perform the following activities:

1. Verify that all products are labeled with the information recorded in the "Label Generation" section. Document the product quantity verified in designated field below.
2. Review the appropriate workstations and surrounding areas for the missing labels.

Investigational Results

☒ N/A

	Label Generation & Reconciliation		Doc No.: FM-0033
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.</i>		Rev.: D
			Page 2 of 2
QC's Initials	GS		Date 07/15/2025



2-WAY 100% SILICONE CLEARTRACT SPT CATHETER

SIZE	16 Fr, 10 mL
QUANTITY	1



Symbols are defined in Instruction for Use (IFU)
 Store away from direct light exposure at room temperature
 Not made with natural rubber latex



SILQ Technologies Corp.
 323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

LOT SLQ-05012025

REF 211610SPT

 2028-05-31



(01) 00850041415154
 (10) SLQ-05012025
 (17) 280531

SLQ-LB-4008 Rev. A

First



2-WAY 100% SILICONE CLEARTRACT SPT CATHETER

SIZE	16 Fr, 10 mL
QUANTITY	1



Symbols are defined in Instruction for Use (IFU)
 Store away from direct light exposure at room temperature
 Not made with natural rubber latex



SILQ Technologies Corp.
 323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

LOT SLQ-05012025

REF 211610SPT

 2028-05-31



(01) 00850041415154
 (10) SLQ-05012025
 (17) 280531

SLQ-LB-4008 Rev. A

Last

NOTES: UNLESS OTHERWISE SPECIFIED

1. GENERAL

1. PRIOR TO PRINTING, VERIFY THE INFORMATION ON THE LABEL IS ACCURATE WHEN COMPARED AGAINST THE WORK ORDER.
2. PRIOR TO PRINTING, VERIFY THE LABEL HAS GONE THROUGH THE INSTALLATION QUALIFICATION (IQ) PROCESS FOR LABELS. REFER TO FM-0044.
3. LABEL CONTENT SHALL BE DEVELOPED PER WI-0009.

2. MATERIAL

1. LABEL: STANDARD WHITE THERMAL TRANSFER
2. RIBBON: WAX / RESIN

3. DIMENSIONS

1. () INDICATES LABEL REFERENCE DIMENSIONS.

4. SAMPLING PLAN

1. IN ADDITION TO THE LABEL GENERATION AND RECONCILIATION FORM, FM-0033, INSPECT N=3 ADDITIONAL SAMPLES AT THE BEGINNING OF THE RUN.

5. GS1 DATA MATRIX VALUES

1. (01): GLOBAL TRADE ITEM NUMBER (GTIN)
2. (10): BATCH OR LOT NUMBER
3. (17): EXPIRATION DATE

6. VISUAL INSPECTION - LABEL TEXT


1. POST-PRINTING, ALL TEXT ON THE LABEL SHALL BE LEGIBLE.
2. ALL VARIABLE LABEL CONTENT (LOT NUMBER, EXPIRATION DATE, ETC.) SHALL MATCH THE WORK ORDER.
3. THE GS1 DATA MATRIX SHALL BE FREE OF MISSING PIXELS AND SHALL NOT HAVE SMUDGES OR STREAKS THROUGH IT. WHEN SCANNED WITH A BARCODE SCANNER, THE TEXT ADJACENT TO THE GS1 DATA MATRIX SHALL MATCH THE TEXT DISPLAYED BY THE SCANNER.
4. (01): VERIFY THE GTIN IS DETAILED IN A HUMAN READABLE FORMAT ADJACENT TO THE GS1 DATA MATRIX. THE GTIN SHALL BE A 14 DIGIT NUMBER.
5. (10): VERIFY THE BATCH OR LOT NUMBER FORMAT IS IN COMPLIANCE WITH QP-019.
6. (17): VERIFY THE EXPIRATION DATE IS IN THE FOLLOWING FORMAT: YYMMDD.
7. VERIFY THE CATHETER STYLE AND CATHETER SIZE ARE REPRESENTATIVE OF A CONFIGURATION FOUND IN THE SKU TABLE.
8. VERIFY THE GTIN IS FOUND IN THE SKU TABLE AND CORRESPONDS WITH THE CORRECT CATHETER SIZE AND STYLE.
9. VERIFY THE REFERENCE NUMBER FOR THE CATHETER CORRESPONDS TO THAT FOUND IN THE SKU TABLE.
10. VERIFY "QUANTITY" IS 1.

7. VISUAL INSPECTION - LABEL SYMBOLS

1. POST-PRINTING, ALL SYMBOLS ON THE LABEL SHALL BE FREE OF SMUDGES / MISSING ARTWORK.
2. THE EXPIRATION DATE NEXT TO ITS CORRESPONDING SYMBOL SHALL BE IN THE FORMAT: YYYY-MM-DD.

REVISIONS			
ZONE	REV.	DESCRIPTION	DATE
	A	INITIAL RELEASE	

ITEM NO.	PARTNO	DESCRIPTION	QTY.
1	SLQ-AW-4008	LABEL, POUCH, SILQ CLEARTRACT SPT CATHETER	1
2	PWY-5007	RIBBON, THERMAL TRANSFER, WAX/RESIN, 4.3W	1
3	SLQ-5004	LABEL, 2.3125 IN. X 7.50 IN.	1



UNLESS OTHERWISE SPECIFIED:
 DIMENSIONS ARE IN MILLIMETERS
 TOLERANCES:
 ANGULAR: MACH ±0.2 BEND ±2.0
 NO PLACE DECIMAL ±0.5
 ONE PLACE DECIMAL ±0.3
 TWO PLACE DECIMAL ±0.15
 THREE PLACE DECIMAL ±0.05

PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PATHWAY NPI. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PATHWAY NPI IS PROHIBITED.



CERTIFIED 13485

DO NOT SCALE DRAWING

pathway
 New Product Introduction


TITLE:
 LABEL, POUCH, SILQ CLEARTRACT SPT CATHETER

SIZE DWG. NO.
B SLQ-LB-4008

REV
 A

SCALE: 1:1

SHEET 1 OF 2



	Label Generation & Reconciliation		Doc No.: FM-0033
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.</i>		Rev.: D
			Page 1 of 2

Label Information					
Part No.	SLQ-LB-4010	Rev.	A	Lot No.	SLQ-05012025 ☐ N/A
Description	Label, Shipper Box, SILQ ClearTract SPT Catheter	UDI No.		Expiration Date:	2028-05-31 ☐ N/A

Label Quantity Calculation	
Required Label Quantity (A)	14
2% Overage (A×0.02=B)	1
QC Retain Labels	2
Quantity of Labels Required (A+B+2=C)	17
Quantity of Labels Printed (C=D)	17
Production's Initial	ST
Date	06/03/2025

Label Verification	
What is the revision of the label per the BOM?	A
What is the revision on the printed label?	A
Is the printed information correct per specifications?	☒ Yes ☐ No
First and last label printed is affixed to back of form	☒ Yes ☐ No
Quantity of labels released to production (E)	15
QC's Initials	AS
Date	06/04/2025

Label Reconciliation	
Quantity of labels used by production (F)	14
Quantity scrapped by production (G)	1
Quantity of unused labels returned to QC and scrapped (H)	0
(E-F-G-H=0)	0
Are all labels accounted for?	☒ Yes ☐ No

If the counts cannot be reconciled, perform the following activities:

1. Verify that all products are labeled with the information recorded in the "Label Generation" section. Document the product quantity verified in designated field below:
2. Review the appropriate workstations and surrounding areas for the missing labels.

Investigational Results

☒ N/A

	Label Generation & Reconciliation		Doc No.: FM-0033
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.</i>		Rev.: D
			Page 2 of 2
QC's Initials	GS	Date	07/15/2025

First



2-WAY 100% SILICONE CLEARTRACT SPT CATHETER

REF

211610SPT

LOT

SLQ-05012025



2028-05-31

QTY

100



SILQ Technologies Corp.

323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

SLQ-LB-4010 Rev. A

Last



2-WAY 100% SILICONE CLEARTRACT SPT CATHETER

REF

211610SPT

LOT

SLQ-05012025



2028-05-31

QTY

100



SILQ Technologies Corp.

323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

SLQ-LB-4010 Rev. A



Label Generation & Reconciliation

Doc No.: FM-0033

This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.

Rev.: D

Page 1 of 2

Label Information

Part No.	SLQ-LB-4010	Rev.	A	Lot No.	SLQ-05012025 ☐ N/A	Label QTY	100 ☐ N/A
Description	Label, Shipper Box, SILQ ClearTract SPT Catheter	UDI No.	☒ N/A			Expiration Date:	2028-05-31 ☐ N/A

Label Quantity Calculation

Required Label Quantity (A)	14		
2% Overage (A×0.02=B)	1		
QC Retain Labels	2		
Quantity of Labels Required (A+B+2=C)	17		
Quantity of Labels Printed (C=D)	17		
Production's Initial	ST	Date	06/12/2025

Label Verification

What is the revision of the label per the BOM?	A		
What is the revision on the printed label?	A		
Is the printed information correct per specifications?	☒ Yes ☐ No		
First and last label printed is affixed to back of form	☒ Yes ☐ No		
Quantity of labels released to production (E)	15		
QC's Initials	GS	Date	06/12/2025

Label Reconciliation

Quantity of labels used by production (F)	14
Quantity scrapped by production (G)	1
Quantity of unused labels returned to QC and scrapped (H)	0
(E-F-G-H=0)	0
Are all labels accounted for?	☒ Yes ☐ No

If the counts cannot be reconciled, perform the following activities:

1. Verify that all products are labeled with the information recorded in the "Label Generation" section. Document the product quantity verified in designated field below.
2. Review the appropriate workstations and surrounding areas for the missing labels.

Investigational Results

☒ N/A

	Label Generation & Reconciliation	Doc No.: FM-0033
	<i>This document contains confidential, proprietary information of Pathway LLC and shall not be copied or reproduced without prior written permission there from.</i>	Rev.: D
		Page 2 of 2
QC's Initials	GS	Date
		07/15/2025



2-WAY 100% SILICONE CLEARTRACT SPT CATHETER

REF 211610SPT

LOT SLQ-05012025

 2028-05-31

QTY 100



SILQ Technologies Corp.

323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

SLQ-LB-4010 Rev. A

First



2-WAY 100% SILICONE CLEARTRACT SPT CATHETER

REF 211610SPT

LOT SLQ-05012025

 2028-05-31

QTY 100



SILQ Technologies Corp.

323 Sunny Isles Blvd, 7th Floor | Sunny Isles Beach, FL 33160

SLQ-LB-4010 Rev. A

Last

NOTES: UNLESS OTHERWISE SPECIFIED

1. GENERAL

1. PRIOR TO PRINTING, VERIFY THE INFORMATION ON THE LABEL IS ACCURATE WHEN COMPARED AGAINST THE WORK ORDER.
2. PRIOR TO PRINTING, VERIFY THE LABEL HAS GONE THROUGH THE INSTALLATION QUALIFICATION (IQ) PROCESS FOR LABELS. REFER TO FM-0044.
3. LABEL CONTENT SHALL BE DEVELOPED PER WI-0009.

2. MATERIAL

1. LABEL: STANDARD WHITE THERMAL TRANSFER
2. RIBBON: WAX / RESIN

3. DIMENSIONS

1. () INDICATES LABEL REFERENCE DIMENSIONS.

4. SAMPLING PLAN

1. IN ADDITION TO THE LABEL GENERATION AND RECONCILIATION FORM, FM-0033, INSPECT N=3 ADDITIONAL SAMPLES AT THE BEGINNING OF THE RUN.

5. VISUAL INSPECTION - LABEL TEXT


1. POST-PRINTING, ALL TEXT ON THE LABEL SHALL BE LEGIBLE.
2. ALL VARIABLE LABEL CONTENT (LOT NUMBER, EXPIRATION DATE, ETC.) SHALL MATCH THE WORK ORDER.
3. VERIFY THE CATHETER STYLE IS REPRESENTATIVE OF A CONFIGURATION FOUND IN THE SKU TABLE.
4. VERIFY THE REFERENCE NUMBER FOR THE CATHETER CORRESPONDS TO THAT FOUND IN THE SKU TABLE.
5. VERIFY "QUANTITY" IS 100.

6. VISUAL INSPECTION - LABEL SYMBOLS

1. POST-PRINTING, ALL SYMBOLS ON THE LABEL SHALL BE FREE OF SMUDGES / MISSING ARTWORK.
2. THE EXPIRATION DATE NEXT TO ITS CORRESPONDING SYMBOL SHALL BE IN THE FORMAT: YYYY-MM-DD.

REVISIONS			
ZONE	REV.	DESCRIPTION	DATE
	A	INITIAL RELEASE	

ITEM NO.	PARTNO	DESCRIPTION	QTY.
1	SLQ-AW-4010	ARTWORK, LABEL, SHIPPER BOX, SILQ CLEARTRACT SPT CATHETER	1
2	PWY-5004	LABEL, 4 X 6	1
3	PWY-5007	RIBBON, THERMAL TRANSFER, WAX/RESIN, 4.3W	1



UNLESS OTHERWISE SPECIFIED:
 DIMENSIONS ARE IN MILLIMETERS
 TOLERANCES:
 ANGULAR: MACH ±0.2 BEND ±2.0
 NO PLACE DECIMAL ±0.5
 ONE PLACE DECIMAL ±0.3
 TWO PLACE DECIMAL ±0.15
 THREE PLACE DECIMAL ±0.05

PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF PATHWAY NPI. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF PATHWAY NPI IS PROHIBITED.



DO NOT SCALE DRAWING



TITLE:
 LABEL, SHIPPER BOX, SILQ CLEARTRACT SPT CATHETER

SIZE DWG. NO.
B SLQ-LB-4010

REV
 A

SCALE: 1:2

SHEET 1 OF 2

PATHWAY INNOVATION LLC

Packing List No: PL-SLQ-0059

Page No: 1 of 1

Company: LGGS INC

Ship Date: 6/16/2025

PO#

Tracking No: 882072221038

Delivery Company: Fedex

Delivery Method: Ground

FOB: Santee, CA

FM-0062, Rev A

STERIS: Ethylene Oxide Certificate Of Processing

Prepared PATHWAY INNOVATION LLC
For
Ethylene Oxide Run ID 1162-18205

<u>Product Code</u>	<u>Lot Number</u>	<u>Quantity (EA)</u>
PWI-7374 Full Cycle	SLQ-05012025	12
Customer ItemID: SLQ-211610SPT		
PWI-7374 Full Cycle	SLQ-05022025	32
Customer ItemID: SLQ-211810SPT		

Preconditioning

<u>Equipment Name</u>	<u>Start Date</u>	<u>End Date</u>	<u>Duration (min)</u>
Preconditioning Room	23-Jun-2025 12:30 AM	23-Jun-2025 5:00 PM	990

		<u>Min</u>	<u>Max</u>
Duration (min)	Spec:	960	1920
Relative Humidity (%RH)	Spec:	40.0	80.0
	Actual:	51.0	63.0
Temperature (°F)	Spec:	105.0	125.0
	Actual:	111.0	118.0

Vessel Processing

<u>Equipment Name</u>	<u>Start Date</u>	<u>End Date</u>	<u>Duration (min)</u>
Chamber 6	23-Jun-2025 5:04 PM	24-Jun-2025 2:02 AM	538
Cycle ID: PWI-7374	Check Value: N/A		

		<u>Min</u>	<u>Max</u>
Duration (min)	Spec:	1	1440

EO Drum Information

<u>Drum ID</u>	<u>Drum Serial Number</u>	<u>Drum Lot Number</u>	<u>EO Usage (LB)</u>
2375	E002727	UTLX930021E25	18.0
Total EO Usage:			18.0

EO Drum Information (Spec)

	<u>EO Usage (LB)</u>
Min	1.0
Max	100.0

Aeration

<u>Equipment Name</u>	<u>Start Date</u>	<u>End Date</u>	<u>Duration (min)</u>	
Aeration Room	24-Jun-2025 2:10 AM	26-Jun-2025 2:15 AM	2885	
		Total Aeration Time:	2885	
			<u>Min</u>	<u>Max</u>
	Duration (min)	Spec:	2880	4320
	Aeration in Room Temperature (°F)	Spec:	103.0	123.0
		Actual:	109.0	115.0

EO Spec ID/Rev/Description ODMS Spec #5 Rev 1 (PWI-7374 Rev. 1 Medical Product)

Product meets Customer specifications; zero nonconformities occurred during this run.

Customer Specification ID: 769

Signature Manifest

Reviewed and E-Signed By: **Antonio Garner (Quality Technician)**

DateTime Esigned 26-Jun-2025 8:03 AM

Document Content Revision 1

Operating facilities are in compliance with applicable state and federal regulations (FDA, NRC, EPA, and OSHA) and provide services under a quality system which meets the requirements of FDA QSR, EN/ISO 13485, and in alignment with EN ANSI/AAMI/ISO 11135.

Processing Location

43425 Business Park Dr
Temecula CA 92590
Phone: 951-694-9340

*** End of Ethylene Oxide Certificate of Processing Document ***

Document ID: 61862

STERIS: Receiving Report
PATHWAY INNOVATION LLC (8130)
Ethylene Oxide

Receipt ID
1161-13286



Product Received From:
PATHWAY INNOVATION LLC
8779 Cottonwood Ave
Suite 105
Santee CA 92071
325

Carrier: NorthStar
Trailer: N/A
Bill of Lading: 3437176
Seal: None
PO(s): P2207

Product counts provided by Customer

Receipt Date/Time: 21-Jun-2025 6:10 PM
Pallet Count: 1
Unloaded By: Brian Neilson
Counted By: Mark Lemus

Container Unload: ☐
Promise Date: 21-Jul-2025 6:10 PM

<u>Item [Pkg]</u>	<u>Lot Number</u>	<u>Quantity</u>			
		<u>Customer</u>	<u>STERIS</u>	<u>CD</u>	<u>UOM</u>
PWI-7374 Full Cycle [1] Customer Item ID: SLQ-211610SPT	SLQ-05012025	12	12	0	EA
PWI-7374 Full Cycle [1] Customer Item ID: SLQ-211810SPT	SLQ-05022025	32	32	0	EA
Receipt Totals:		44	44	0	

Inventory by Received Pallet

<u>Pallet (SN)</u>	<u>Item ID [Pkg]</u>	<u>Lot Number</u>	<u>Quantity</u>			<u>Warehouse</u>	<u>Pallet</u>
			<u>Customer</u>	<u>STERIS</u>	<u>UOM</u>	<u>Location</u>	<u>Total</u>
1 (103951)	PWI-7374 Full Cycle [1] Customer Item ID: SLQ-211610SPT	SLQ-05012025	12	12	EA	N/A	1
1 (103951)	PWI-7374 Full Cycle [1] Customer Item ID: SLQ-211810SPT	SLQ-05022025	32	32	EA	N/A	1

Document ID: 61714

STERIS: Receiving Report
PATHWAY INNOVATION LLC (8130)
Ethylene Oxide

Receipt ID
1161-13286



Product Received From:
PATHWAY INNOVATION LLC
8779 Cottonwood Ave
Suite 105
Santee CA 92071
325

Carrier: **NorthStar**
Trailer: **N/A**
Bill of Lading: **3437176**
Seal: **None**
PO(s): **P2207**

Product counts provided by Customer

Receipt Date/Time: **21-Jun-2025 6:10 PM**

Pallet Count: **1**

Unloaded By: **Brian Neilson**

Counted By: **Mark Lemus**

Container Unload:



Promise Date: **21-Jul-2025 6:10 PM**

Signature Manifest

Reviewed By: **Mark Lemus (Material Handler II)**

Date/Time Signed 21-Jun-2025 6:13 PM

Document Content Revision 1

Processing Location

43425 Business Park Dr
Temecula CA 92590
Phone: 951-694-9340

Document ID: 61714

PL 6-24-25


QA Reviewed

AG JUN 24 2025

HONEYWELL
30755312-001

Tulane
FT. WASHINGTON, P.A. U.S.

6/23/25





Applied Sterilization Technologies

EO Process Report

STERIS

Applied Sterilization Technologies

43425 Business Park Drive

Temecula, CA 92590

Tel: (951) 694-9340

Fax: (951) 699-5981

Customer Name Pathway FC

Load # PWI-1162-18205

Date Processed 06/23/2025 17:04:31

Pressure Units "HgA

Cycle Number 7374

Cycle Last Modified: 4/30/2018 1:56:49PM

Chamber # 6

Cycle Counter # 5,122

Temperature Units Fahrenheit

Phase Name

Initial Evacuation

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/23/2025 17:04:31	00:00:00	29.0	127
06/23/2025 17:09:31	00:05:00	21.6	125
06/23/2025 17:14:32	00:10:00	14.8	125
06/23/2025 17:19:32	00:15:00	8.1	124
06/23/2025 17:24:32	00:20:00	3.0	124
06/23/2025 17:26:50	00:22:19	2.0	124

Summary of Initial Evacuation

Pressure Start	29.0	Min Pressure	2.0	Min Temp	124	Pressure Change	-27.0
Pressure End	2.0	Max Pressure	29.0	Max Temp	127	Phase Time	00:22:19

Stabilization

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/23/2025 17:26:50	00:22:19	2.0	124
06/23/2025 17:28:50	00:24:19	2.1	125
06/23/2025 17:28:50	00:24:19	2.1	125

Summary of Stabilization

Pressure Start	2.0	Min Pressure	2.0	Min Temp	124	Pressure Change	0.1
Pressure End	2.1	Max Pressure	2.1	Max Temp	125	Phase Time	00:02:00

Leak Check

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/23/2025 17:28:50	00:24:19	2.1	125
06/23/2025 17:29:51	00:25:19	2.1	126
06/23/2025 17:30:50	00:26:19	2.2	126
06/23/2025 17:31:51	00:27:19	2.2	126
06/23/2025 17:32:51	00:28:19	2.2	126
06/23/2025 17:33:50	00:29:19	2.2	126

Summary of Leak Check

Pressure Start	2.1	Min Pressure	2.1	Min Temp	125	Pressure Change	0.1
Pressure End	2.2	Max Pressure	2.2	Max Temp	126	Phase Time	00:05:00

Initial Nitrogen Wash

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/23/2025 17:33:50	00:29:19	2.2	126
06/23/2025 17:38:51	00:34:19	9.3	128
06/23/2025 17:43:51	00:39:19	16.8	128
06/23/2025 17:48:51	00:44:19	24.3	128
06/23/2025 17:51:21	00:46:49	28.0	128

Customer Name Pathway FC
Load # PWI-1162-18205

STERIS
Applied Sterilization Technologies
Temecula, CA

Summary of Initial Nitrogen Wash

Pressure Start	2.2	Min Pressure	2.2	Min Temp	126	Pressure Change	25.8
Pressure End	28.0	Max Pressure	28.0	Max Temp	128	Phase Time	00:17:30

Initial Evacuation

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/23/2025 17:51:21	00:46:49	28.0	128
06/23/2025 17:56:20	00:51:49	21.3	124
06/23/2025 18:01:21	00:56:50	14.3	123
06/23/2025 18:06:21	01:01:50	7.3	123
06/23/2025 18:11:21	01:06:50	2.8	123
06/23/2025 18:13:08	01:08:37	2.0	123

Summary of Initial Evacuation

Pressure Start	28.0	Min Pressure	2.0	Min Temp	122	Pressure Change	-26.0
Pressure End	2.0	Max Pressure	28.0	Max Temp	128	Phase Time	00:21:47

Humidification

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp	Avg. % RH
06/23/2025 18:13:08	01:08:37	2.0	123	0
06/23/2025 18:14:09	01:09:37	2.6	124	12
06/23/2025 18:15:08	01:10:37	2.9	125	27
06/23/2025 18:15:19	01:10:48	3.0	125	28

Summary of Humidification

Pressure Start	2.0	Min Pressure	2.0	Min Temp	123	Pressure Change	1.0	Min RH	0
Pressure End	3.0	Max Pressure	3.0	Max Temp	125	Phase Time	00:02:11	Max RH	28

Humidity Dwell

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp	Avg. % RH
06/23/2025 18:15:19	01:10:48	3.0	125	28
06/23/2025 18:20:19	01:15:48	3.0	125	32
06/23/2025 18:25:19	01:20:48	3.0	125	32
06/23/2025 18:30:19	01:25:48	3.0	125	31
06/23/2025 18:35:20	01:30:48	3.0	125	31
06/23/2025 18:40:20	01:35:48	3.0	125	30
06/23/2025 18:45:20	01:40:48	3.0	125	30
06/23/2025 18:50:19	01:45:48	3.0	125	30
06/23/2025 18:55:20	01:50:48	3.0	125	30
06/23/2025 19:00:20	01:55:49	3.1	125	30
06/23/2025 19:05:20	02:00:49	3.1	125	30
06/23/2025 19:10:20	02:05:49	3.1	125	30
06/23/2025 19:15:19	02:10:48	3.1	125	31

Summary of Humidity Dwell

Pressure Start	3.0	Min Pressure	3.0	Min Temp	125	Pressure Change	0.1	Min RH	28
Pressure End	3.1	Max Pressure	3.1	Max Temp	125	Phase Time	01:00:00	Max RH	32

Sterilant Injection

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp	EO Drum 1 (lb)	EO Drum 2 (lb)	EO CONC
06/23/2025 19:15:19	02:10:48	3.1	125	537	0	0
06/23/2025 19:20:20	02:15:48	7.7	126	532	0	253
06/23/2025 19:25:20	02:20:48	12.8	126	524	0	563
06/23/2025 19:26:55	02:22:24	14.3	126	522	0	649

Customer Name Pathway FC
Load # PWI-1162-18205

STERIS
Applied Sterilization Technologies
Temecula, CA

Summary of Sterilant Injection

Pressure Start 3.1	Min. Pressure 3.1	Min. Temp 125	Drum 1 Start 537	Drum 2 Start 0
Pressure End 14.3	Max Pressure 14.3	Max Temp 126	Drum 1 End 522	Drum 2 End 0
			EtO Used 1 15	EtO Used 2 0
Pressure Change 11.2	MIN EO CONC 0.00			
Phase Time 00:11:36	MAX EO CONC 649.00		Total EO Weight 15	

Exposure

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp	EO Drum 1 (lb)	EO Drum 2 (lb)	EO CONC
06/23/2025 19:26:55	02:22:24	14.3	126	522	0	649
06/23/2025 19:31:55	02:27:24	14.3	125	523	0	673
06/23/2025 19:36:55	02:32:24	14.2	125	522	0	661
06/23/2025 19:41:56	02:37:24	14.3	125	523	0	665
06/23/2025 19:46:56	02:42:24	14.2	125	522	0	656
06/23/2025 19:51:56	02:47:24	14.2	125	522	0	649
06/23/2025 19:56:56	02:52:25	14.2	125	522	0	643
06/23/2025 20:01:56	02:57:25	14.1	125	522	0	638
06/23/2025 20:06:56	03:02:25	14.1	125	521	0	634
06/23/2025 20:11:56	03:07:25	14.3	125	522	0	642
06/23/2025 20:16:56	03:12:25	14.3	125	521	0	639
06/23/2025 20:21:57	03:17:25	14.3	125	521	0	637
06/23/2025 20:26:57	03:22:25	14.3	125	521	0	633
06/23/2025 20:31:57	03:27:25	14.3	125	521	0	630
06/23/2025 20:36:56	03:32:25	14.3	125	521	0	627
06/23/2025 20:41:56	03:37:25	14.3	125	520	0	624
06/23/2025 20:46:57	03:42:25	14.3	125	520	0	623
06/23/2025 20:51:57	03:47:25	14.3	125	520	0	621
06/23/2025 20:56:57	03:52:26	14.3	125	520	0	619
06/23/2025 21:01:57	03:57:26	14.3	125	520	0	619
06/23/2025 21:06:57	04:02:26	14.3	125	520	0	617
06/23/2025 21:11:57	04:07:26	14.3	125	520	0	614
06/23/2025 21:16:57	04:12:26	14.2	125	520	0	615
06/23/2025 21:21:57	04:17:26	14.2	125	520	0	613
06/23/2025 21:26:57	04:22:26	14.2	125	520	0	613
06/23/2025 21:31:58	04:27:26	14.2	125	520	0	613
06/23/2025 21:36:58	04:32:26	14.2	125	520	0	611
06/23/2025 21:41:57	04:37:26	14.2	125	520	0	611
06/23/2025 21:46:57	04:42:26	14.2	125	520	0	611
06/23/2025 21:51:58	04:47:26	14.2	125	520	0	608
06/23/2025 21:56:58	04:52:26	14.2	125	519	0	608
06/23/2025 22:01:58	04:57:27	14.2	125	519	0	608
06/23/2025 22:06:58	05:02:27	14.2	125	519	0	606
06/23/2025 22:11:59	05:07:27	14.2	125	519	0	607
06/23/2025 22:16:58	05:12:27	14.2	125	519	0	607
06/23/2025 22:21:58	05:17:27	14.2	125	519	0	604
06/23/2025 22:26:58	05:22:27	14.2	125	519	0	604
06/23/2025 22:31:58	05:27:27	14.2	125	519	0	604
06/23/2025 22:36:59	05:32:27	14.2	125	519	0	602
06/23/2025 22:41:59	05:37:27	14.2	125	519	0	602
06/23/2025 22:46:59	05:42:27	14.2	125	519	0	601
06/23/2025 22:51:58	05:47:27	14.2	125	519	0	601
06/23/2025 22:56:59	05:52:27	14.2	125	519	0	601
06/23/2025 23:01:59	05:57:28	14.2	125	519	0	601
06/23/2025 23:06:59	06:02:28	14.2	125	519	0	600
06/23/2025 23:11:59	06:07:28	14.2	125	519	0	598
06/23/2025 23:16:59	06:12:28	14.2	125	519	0	600
06/23/2025 23:21:59	06:17:28	14.2	125	519	0	599
06/23/2025 23:26:56	06:22:24	14.2	125	519	0	597

Customer Name Pathway FC
Load # PWI-1162-18205

STERIS
 Applied Sterilization Technologies
 Temecula, CA

Summary of Exposure

Pressure Start 14.3	Min Pressure 14.1	Min Temp 125	Drum 1 Start 522	Drum 2 Start 0
Pressure End 14.2	Max Pressure 14.4	Max Temp 126	Drum 1 End 519	Drum 2 End 0
			EtO Used 1 3	EtO Used 2 0
Pressure Change -0.1		MIN EO CONC 597.00		
Phase Time 04:00:00		MAX EO CONC 675.00		
			Total EO Weight 3	

Sterilant Removal

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/23/2025 23:26:56	06:22:24	14.2	125
06/23/2025 23:31:56	06:27:24	9.9	124
06/23/2025 23:36:56	06:32:24	5.2	123
06/23/2025 23:41:55	06:37:24	2.2	123
06/23/2025 23:42:33	06:38:02	2.0	123

Summary of Sterilant Removal

Pressure Start 14.2	Min Pressure 2.0	Min Temp 123	Pressure Change -12.2
Pressure End 2.0	Max Pressure 14.2	Max Temp 125	Phase Time 00:15:38

Post Nitrogen Wash

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/23/2025 23:42:33	06:38:02	2.0	123
06/23/2025 23:47:33	06:43:02	9.8	126
06/23/2025 23:52:34	06:48:02	17.0	127
06/23/2025 23:57:34	06:53:02	24.3	127
06/24/2025 00:00:06	06:55:34	28.0	127

Summary of Post Nitrogen Wash

Pressure Start 2.0	Min Pressure 2.0	Min Temp 123	Pressure Change 26.0
Pressure End 28.0	Max Pressure 28.0	Max Temp 127	Phase Time 00:17:32

Post Evacuation

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/24/2025 00:00:06	06:55:34	28.0	127
06/24/2025 00:05:06	07:00:34	21.3	123
06/24/2025 00:10:05	07:05:35	14.3	122
06/24/2025 00:15:06	07:10:35	7.4	122
06/24/2025 00:20:06	07:15:35	2.9	122
06/24/2025 00:22:19	07:17:48	2.0	122

Summary of Post Evacuation

Pressure Start 28.0	Min Pressure 2.0	Min Temp 122	Pressure Change -26.0
Pressure End 2.0	Max Pressure 28.0	Max Temp 127	Phase Time 00:22:13

Post Nitrogen Wash

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/24/2025 00:22:19	07:17:48	2.0	122
06/24/2025 00:27:19	07:22:48	9.1	126
06/24/2025 00:32:19	07:27:48	16.6	127
06/24/2025 00:37:20	07:32:48	24.1	127
06/24/2025 00:39:59	07:35:28	28.0	127

Summary of Post Nitrogen Wash

Pressure Start 2.0	Min Pressure 2.0	Min Temp 122	Pressure Change 26.0
Pressure End 28.0	Max Pressure 28.0	Max Temp 128	Phase Time 00:17:40

Customer Name Pathway FC
Load # PWI-1162-18205

STERIS
Applied Sterilization Technologies
Temecula, CA

Post Evacuation

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/24/2025 00:39:59	07:35:28	28.0	127
06/24/2025 00:45:00	07:40:28	21.3	123
06/24/2025 00:50:00	07:45:28	14.3	122
06/24/2025 00:55:00	07:50:28	7.3	122
06/24/2025 01:00:00	07:55:28	2.8	122
06/24/2025 01:01:55	07:57:23	2.0	122

Summary of Post Evacuation

Pressure Start	28.0	Min Pressure	2.0	Min Temp	121	Pressure Change	-26.0
Pressure End	2.0	Max Pressure	28.0	Max Temp	127	Phase Time	00:21:55

Post Nitrogen Wash

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/24/2025 01:01:55	07:57:23	2.0	122
06/24/2025 01:06:55	08:02:23	9.1	126
06/24/2025 01:11:55	08:07:23	16.6	127
06/24/2025 01:16:54	08:12:24	24.1	127
06/24/2025 01:19:34	08:15:03	28.0	128

Summary of Post Nitrogen Wash

Pressure Start	2.0	Min Pressure	2.0	Min Temp	122	Pressure Change	26.0
Pressure End	28.0	Max Pressure	28.0	Max Temp	128	Phase Time	00:17:40

Post Evacuation

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/24/2025 01:19:34	08:15:03	28.0	128
06/24/2025 01:24:34	08:20:03	21.3	124
06/24/2025 01:29:35	08:25:03	14.3	122
06/24/2025 01:34:35	08:30:03	7.3	122
06/24/2025 01:39:35	08:35:03	2.8	122
06/24/2025 01:41:22	08:36:50	2.0	122

Summary of Post Evacuation

Pressure Start	28.0	Min Pressure	2.0	Min Temp	121	Pressure Change	-26.0
Pressure End	2.0	Max Pressure	28.0	Max Temp	128	Phase Time	00:21:47

Air Backfill

Date / Time	Cycle Time	Avg. CH Press	Avg. CH Temp
06/24/2025 01:41:22	08:36:50	2.0	122
06/24/2025 01:46:21	08:41:50	7.8	125
06/24/2025 01:51:22	08:46:51	14.1	125
06/24/2025 01:56:22	08:51:51	20.4	125
06/24/2025 02:01:22	08:56:51	26.5	124
06/24/2025 02:02:35	08:58:04	27.9	124

Summary of Air Backfill

Pressure Start	2.0	Min Pressure	2.0	Min Temp	122	Pressure Change	25.9
Pressure End	27.9	Max Pressure	27.9	Max Temp	126	Phase Time	00:21:13

Customer Name Pathway FC
Load # PWI-1162-18205

STERIS
Applied Sterilization Technologies
Temecula, CA

End of Cycle

Production Review

Ruby Rogers
Signature

Date


6-24-25

Quality Review

[Signature]
Signature

Date

6-26-25



QA Reviewed
JUN 28 2025

pn 6-24-25

pen G-25-25

Date Inspected	Part No.	Rev.	Description
07/15/2025	SLQ-211610SPT	B	SILQ 2-way Cleartract SPT 16 Fr x 10 ml Catheter, Sterile, 10-Pack
Lot No.	Lot QTY		
SLQ-05012025	1298 catheters/ 129.8 shippers		
Specification	Identification of the product is traceable between the DHR, sterilization records, and laboratory report	DHR is completed and meets good documentation practices per QP-003	Sterilization processing record meets parameters of validated process
QP-010 Revision	F	F	F
Lower Spec Limit (LSL)	N/A	N/A	N/A
Upper Spec Limit (USL)	N/A	N/A	N/A
UOM	N/A	N/A	N/A
Measuring Equipment Type	Visual	Visual	Visual
Asset No.	N/A	N/A	N/A
Cal Due Date	N/A	N/A	N/A
AQL	See Below	See Below	See Below
Sample Size	DHR Sterilization Processing Record Laboratory Report	DHR	Sterilization Processing Record
Sample No.	Variable/ Attribute Data		
1	1	1	1
Pass QTY	1	1	1
Fail QTY	0	0	0
Result (P/F)	P	P	P
Min	N/A	N/A	N/A
Max	N/A	N/A	N/A
Mean	N/A	N/A	N/A
StDev	N/A	N/A	N/A
Comments	None		
Disposition:	Accept	NCR No:	None
Inspected By:	Larry Dacorn	Date:	07/15/2025
Verified By:	Larry Dacorn	Date:	07/15/2025