


ClearTract™

FOLEY CATHETER



FDA CLEARED

For Suprapubic, Urethral, and
Nephrostomy Use



PUBLISHED MAY 2022 IN ADVANCED MATERIALS

SCAN
HERE
TO
READ



Silq's groundbreaking zwitterion surface technology

For the past twenty years, scientists have shown that bio-inspired polymers known as polyzwitterions strongly resist the deposition of microbes and proteins on surfaces that lead to biofilm formation.⁵⁻¹⁵ Silq's proprietary polyzwitterionic surface technology was designed to reduce the risk of infection and other complications associated with implanted medical devices.¹⁶ Read more about this in our *Advanced Materials* article (left).

Silq has been able to successfully apply zwitterionic polymers to medical devices at commercial scale, finally unlocking the tremendous potential of zwitterions. Additionally, Silq's technology does not rely on the use of antibiotic drugs or silver ions that are known to kill microbes in the body. These antimicrobials can lead to patient discomfort¹⁷ and the dangerous proliferation of "superbugs" that are resistant to modern antibiotic treatments.¹⁸⁻²⁰

ClearTract™ Catheters: Improves Patient Comfort and Safety

SILQ'S TREATMENT INCREASES CATHETER LUBRICITY BY UP TO 18X

MADE FROM
100% MEDICAL
GRADE SILICONE

NO
LATEX

NO
BPA


NO
DEHP

The days of latex healthcare products are coming to an end - 4.3% of the world population has a sensitivity to latex. Healthcare workers have more than double the latex sensitivity rate due to repeated exposure.²⁴

ClearTract™ catheters are made from 100% medical grade silicone, protecting both patients and their healthcare providers.

Regarding safety, ClearTract™ catheters have been cleared by the FDA based on the most recent and more rigorous biocompatibility testing requirements for catheters.

DRAG FORCE (GRAMS)



TESTING SUCCESSFULLY COMPLETED

- Cytotoxicity
- Systemic Toxicity
- Irritation
- Pyrogenicity
- Sub-acute Toxicity
- Sensitivity
- Sub-chronic Toxicity
- Genotoxicity*

*Additional requirement specifically for suprapubic clearance

**All test data available on file

ClearTract™ Foley Catheter - Ordering Information

CATHETER SIZE	BALLOON SIZE	VALVE	CATHETER SKU#
14F	10cc	2-Way	21400101003
16F	10cc	2-Way	21600101003
18F	10cc	2-Way	21800101003

REFERENCES:

1) B Foxman. **The epidemiology of urinary tract infection.** Nat Rev Urol. 2010, 7(12): 653–660. 2) C S Hollenbeak, et al. **The attributable cost of catheter-associated urinary tract infections in the United States: A systematic review.** Am J Infect Control. 2018, 46(7): 751-757. 3) CDC. **Antibiotic Resistance Threats in the United States, 2019.** Atlanta, GA: U.S. Department of Health and Human Services, CDC; 2019. 4) Lo, Evelyn, et al. **"Strategies to Prevent Catheter-Associated Urinary Tract Infections in Acute Care Hospitals: 2014 Update."** Infection Control & Hospital Epidemiology, vol. 35, no. 5, May 2014, pp. 464–479. 5) Singha, Priyadarshini, et al. **"A Review of the Recent Advances in Antimicrobial Coatings for Urinary Catheters."** Acta Biomaterialia, vol. 50, Mar. 2017, pp. 20–40. 6) Schleno JB. **Zwitterion: coating surfaces with zwitterionic functionality to reduce nonspecific adsorption.** Langmuir. 2014 Aug 19;30(32):9625-36. 7) Colilla M, **The Role of Zwitterionic Materials in the Fight against Proteins and Bacteria.** Medicines. 2018 Nov 22;5(4):125. 8) Kazuhiko Ishihara. **Biomimetic materials based on zwitterionic polymers toward human-friendly medical devices,** Science and Technology of Advanced Materials, 23:1, 498-524. 9) Liu S **Recent Advances in Zwitterionic Hydrogels: Preparation, Property, and Biomedical Application.** Gels. 2022; 8(1):46. 10) Yanxian Zhang et al. **Fundamentals and applications of zwitterionic antifouling polymers.** 2019 J. Phys. D: Appl. Phys. 52 403001. 11) Li Q, **Zwitterionic Biomaterials.** Chem Rev. 2022 Dec 14;122(23):17073-17154. 12) Sin, MC. **Hemocompatibility of zwitterionic interfaces and membranes.** Polym J 46, 436–443 (2014). 13) Racovita S. **Polybetaines in Biomedical Applications.** International Journal of Molecular Sciences. 2021; 22(17):9321. 14) Ladd, Jon. **Zwitterionic Polymers Exhibiting High Resistance to Nonspecific Protein Adsorption from Human Serum and Plasma.** Biomacromolecules (2008). 9. 1357-61. 15) Zhang, L, Cao, Z, Bai, T. et al. **Zwitterionic hydrogels implanted in mice resist the foreign-body reaction.** Nat Biotechnol 31, 553–556 (2013). 16) B McVerry, et al. **A Readily Scalable, Clinically Demonstrated, Antibiofouling Zwitterionic Surface Treatment for Implantable Medical Devices.** Adv. Mater. 2022, 34(20): 2200254. 17) Pickard R. **Types of urethral catheter for reducing symptomatic urinary tract infections in hospitalised adults requiring short-term catheterisation: multicentre randomized controlled trial and economic evaluation of antimicrobial- and antiseptic impregnated urethral catheters.** Health Technol Assess 2012;16(47). 18) McNeilly, Oliver, et al. **"Emerging Concern for Silver Nanoparticle Resistance in Acinetobacter Baumannii and Other Bacteria."** Frontiers in Microbiology, vol. 12, 16 Apr. 2021. 19) Hosny, Alaa El-Dien MS, et al. **"The Increasing Threat of Silver-Resistance in Clinical Isolates from Wounds and Burns."** Infection and Drug Resistance, vol. Volume 12, July 2019, pp. 1985–2001. 20) Köves B, Magyar A, Tenke P. **Spectrum and antibiotic resistance of catheter-associated urinary tract infections.** GMS Infect Dis. 2017 06.

SILQ Technologies Corp.™

www.silq.tech

info@silq.tech

(424) 309-8523

323 Sunny Isles Blvd, 7th Floor Sunny Isles Beach, FL 33160

DISTRIBUTED BY
FRONT1
HEALTH PARTNERS

info@front1hp.com (866) 277-0124