

Silq Technologies Corp. Management Review

Timeframe under review: 2024



Meeting date & attendees

- Meeting date: 11 Feb 2025
- Meeting start time: 1 PM PST
- Attendees and roles:
 - Verne Sharma, President & CEO
 - Brian McVerry, Chief Technology Officer
 - Tom Downey, Chief Financial Officer
 - Chuck Greiner, Chief Commercial Officer
 - Sigi Caron, VP Regulatory & Clinical Affairs
 - Tavi Cruz-Uribe, VP Quality Assurance (Management Representative)
 - Chris Turner, Director, Manufacturing
 - Bryan Wong, Director, QA
 - Ethan Rao, Director, R&D
 - Diana Porter, Quality Specialist (Independent reviewer)



Agenda

- Review of minutes and action items from previous management review meeting
- Review of the Quality Policy
- Review of Quality Objectives
- Review of audits (Internal audits, external audits, and supplier audits)
- Status of the Corrective and Preventive Action System
- Review of product complaints trends
- Review of production data including nonconformities
- Review of training system performance
- Review of Post Market Surveillance / customer feedback data
- Review of reportings to regulatory authorities
- Review of statistical techniques
- Review of supplier performance
- Changes that could affect the quality management system
- Recommendations for improvement
- New or revised standards and regulatory requirements
- Review of Risk Management Activities
- Review status of Design Control Activities
- Review of Resource Requirements



Review of minutes and action items from previous Management Review meeting



Review of minutes and action items from previous MR meeting

- Review of 2023 Management Review meeting minutes (held on 03 Feb 2024)
- Review of action items from 2023 MR meeting: 10 action items from 2023. 6 action items have been completed. 2 action items were addressed through a different corrective action. 2 action items postponed.

#	Topic	Closure details	Status
1	Review of minutes and action items from the previous meeting	Action item: Complete action items 2, 4, and 5 from 2022 Management Review meeting. Quality objectives established in QM.SLQ037 Rev A. Template for complaint response to complainant released (TMP1-QMSLQ021). Quality agreement with Pathway updated and signed on 04 Mar 2024.	Completed
2	Review of the Quality Policy	Action item: Update quality policy to add "Excellent" before "Patient care..." Quality policy updated in QM.SLQ027 Rev D and QM.SLQ035 Rev C.	Completed
3	Review of general quality goals and objectives	Action item: Release quality objectives. Quality objectives released in QM.SLQ037 Rev A.	Completed



Review of minutes and action items from previous MR meeting

#	Topic	Closure details	Status
4	Review of general quality goals and objectives	Action item: Assess complaint rate of Silq catheters from data obtained in 2022-2023 to determine if we meet 5% complaint rate quality objective. Differentiate from reported device deficiencies in the study Based on an assessment performed in Feb 2024, the complaint rate is 0.55%	Completed
5	Review of Product Complaint trends	Action item: Update catheter specifications to include gap specification of 1mm or less. Ningbo implemented a new valve that ensures placement and prevents dislodging.	Addressed through a design change
6	Review of Product Complaint trends	Action item: Add incoming inspection requirement to measure the gap. Ningbo implemented a new valve that ensures placement and prevents dislodging.	Addressed through a design change
7	Status of the Production Data including nonconformities	Action item: Confirm with Pathway the shelf life for all product lines is 3 years. Pathway to update labels, accordingly. Confirmed with Pathway all product lines are at a 3 year shelf life. All finished goods inventory have been re-labeled accordingly.	Completed



Review of minutes and action items from previous MR meeting

#	Topic	Closure details	Status
8	General recommendations for improvement	Action item: Complete Gage R&R study on in-process quality control test methods for UV Light Absorbance and Rinse Test.	Postponed
9	General recommendations for improvement	Action item: Complete test method validation for 1H and 19F NMR spectra readings.	Postponed
10	General recommendations for improvement	Action item: Extend shelf-life of suspension. Suspension shelf-life has been extended to two (2) years. Reference report VV.SLQ029.	Completed



Review of Quality Policy



Review of Quality Policy

Excellent patient care is the focus of all SILQ activities. We endeavor to meet or exceed our customers' satisfaction and expectations in providing safe and effective coated medical devices of the highest quality and that meets all applicable regulatory requirements.

Each SILQ employee is personally committed to achieving this goal and continually improving our ability to maintain customer focus by operating a highly effective quality system that meets the requirements of applicable global standards and regulatory requirements for medical devices.



Review of Quality Objectives



Quality Objectives

1. Objective: Incoming quality

- a. **Goal:**
The lot acceptance rate of incoming material used in finished product shall be greater than 80%.
- b. **Results: Objective met**
Acceptance rate = 100%. Two (2) lots received in 2024. Both lots accepted.
 - SLQ-6000 Catheter, 14 Fr x 10 mL, Silicone, NY Foley, 2-Way
 - SLQ-6002 Catheter, 18 Fr x 10 mL, Silicone, NY Foley, 2-Way

2. Objective: Finished product quality

- a. **Goal:**
The complaint rate of failed catheters shall be less than 5% of the overall distributed catheters.
- b. **Results: Objective met**
The complaint rate in 2024 was 0.3%. A total of 1,737 catheters (from all sizes) were sold in 2024. Of the four customer complaints received in 2024, there were a total of at least five failed catheters.



Review of Audits

(Internal, External, Supplier)



Summary of audits

Internal audit

Audit	Scope of Audit	Date(s)	# Findings/AOC	Corrective Actions	Status of CA's
Internal (#24-01)	ISO 13485:2016, 21 CFR Parts 803, 806 & 820	January 16-17 th , 2025	0 findings 0 recommendations 0 concerns	N/A	N/A



Summary of audits (cont.)

Upcoming 2025 Audits

1. Internal audit schedule 2025

Internal Quality Audit Schedule		Year: <u>2025</u>				
Title	Standard / Regulation	Q1	Q2	Q3	Q4	Date Completed
Quality Management System	ISO 13485 Section 4				X	
Management Responsibility	ISO 13485 Section 5				X	
Resource Management	ISO 13485 Section 6				X	
Product Realization	ISO 13485 Section 7				X	
Measurement, Analysis, and Improvement	ISO 13485 Section 8				X	
Medical Device Reporting (FDA)	21 CFR Part 803				X	
Reports of Corrections and Removals (FDA)	21 CFR Part 806				X	

"X" indicates when audit is scheduled to occur. Audits must take place within the designated quarter.



Status of CAPA System



CAPA System Summary

0 CAPA's initiated in 2024



Review of Product Complaints



Review of Product Complaints

4 product complaints were received in 2024



Review of Product Complaints

PCF# 2024001

- **Complaint status: Closed**
- Complainant: Oliver Estrada, Ranchos Los Amigos National Rehabilitation Center
- Location: Downey, CA
- Date reported: March 11, 2024
- Number of units: > 2 (16Fr x 10mL)
- Complaint description: Two (2) ClearTract catheters had a dark spot on the distal end of the catheter itself and is visible through the clear sterile packaging.
- Patient injury: None
- Reportability determination: Not a reportable complaint
- Failure investigation: Insufficient cleaning of the machinery which led to a build-up of Silq's surface treatment polymer in the underlying manifold which resulted in surface agglomerates.
- Corrective action: Pathway implemented a visual inspection of coated catheters to catch any with surface defects and agglomerations.



Review of Product Complaints

PCF# 2024002

- **Complaint status: Closed**
- Complainant: Tousif Hussain Sheikh, HPFY
- Location: N/A
- Date reported: April 3, 2024
- Number of units: 1 (18Fr x 10mL)
- Complaint description: the ClearTract 18 Fr x 10ml catheter (Model #211810SPT) was defective with the issue that "she recently went to the doctor and right when they were trying to pump the little balloon, one of the things where the syringe goes in broke."
- Patient injury: None
- Reportability determination: Not a reportable complaint
- Failure investigation: Defective catheter was not returned for investigation



Review of Product Complaints

PCF# 2024003

- **Complaint status: Closed**
- Complainant: Gina Koehler, Cleveland Clinic
- Location: Cleveland, OH
- Date reported: May 22, 2024
- Number of units: 1 (16Fr x 10mL)
- Complaint description: When the nurse was filling the balloon with saline, the plastic end of the tube where the syringe is inserted came off. As a result, she was unable to adequately fill the balloon and we had to use another catheter..
- Patient injury: None
- Reportability determination: Not a reportable complaint
- Failure investigation: The catheter was returned to Silq for analysis, however the Luer lock connector which malfunctioned was not present as it had been discarded at the time of the malfunction. After discussion with the complainant, when the syringe containing saline was connected to the Luer lock connector, the Luer lock became unseated from the catheter neck.
- Corrective action: Ningbo implemented a new valve to ensure proper placement and prevents dislodging.



Review of Product Complaints

PCF# 2023004

- **Complaint status: Closed**
- Complainant: Gurav Pandey, HPFY
- Location: N/A
- Date reported: December 17, 2024
- Number of units: 1 (18Fr x 10mL)
- Complaint description: According to the customer, one catheter (item #21800101003, lot #SLQ-81020515241) was defective and had a crack.
- Patient injury: None
- Reportability determination: Not a reportable complaint
- Failure investigation: Aside from normally observed encrustation on Foley catheters, no other apparent defect was observed. No crack or tear could be found on the catheter after a thorough evaluation of the catheter. There was also no apparent issue with the tip or valve of the catheter based on a visual inspection. Additional testing was performed to confirm if there was a crack or leak in the device. The catheter balloon was filled with 5mL of water for 1 minute before being increased to 10mL for an additional 1 minute. No issues or leaks were observed.



Review of Production Data



Review of Production Data

2024 Build Summary

Finished devices

Assembly	Lot #	Qty released	Exp Date
SLQ-211810SPT	SLQ-11182024	0 (lot failed LAL, 60.3 cases)	30 Nov 2027
SLQ-211410SPT	SLQ-11192024	22.1 cases	30 Nov 2027
SLQ-211610SPT	SLQ-11202024	94.4 cases	30 Nov 2027
SLQ-211810SPT	SLQ-81020515241	38 cases	15 May 2027

Suspension

Assembly	Lot #	Qty released	Mfg Date	Exp Date
C.SLQ001	None	--	--	--



Review of Production Data (cont.)

Nonconformance Summary 2024

- No NCR's were initiated by Silq or Pathway in 2024



Review of Training System performance



Review of Training System

- Training Status (based on employee's training plan):
 - All employee training is up to date.
 - There were no new hires in 2024.



Review of Post Market Surveillance / Customer Feedback Data



Review of Post Market Surveillance / Customer Feedback Data

Post-Market Surveillance Data	Summary of activity/results
Post-market Clinical Trials	Analysis report on the primary endpoint of the RCT Study (analysis of biofilm) was submitted by Technomics on 19 Dec 2024. Results indicate that the ClearTract had significantly less biofilm than the BARDEX and Rusch Gold Standard catheters using both the first and second models.



Review of Post Market Surveillance / Customer Feedback Data

Post-Market Surveillance Data	Summary of activity/results
MAUDE database	Database search from Jan 1, 2024 to Dec 31, 2024 for reports in the following Product Classes. • Catheter, Urological: ○ 1,275 reports ▪ 1,135 malfunctions, 139 patient injuries, 1 death • Catheter, Urological (Antimicrobial) and Accessories ○ 1,978 reports ▪ 1,914 malfunctions, 64 patient injuries, 0 deaths <u>Common malfunctions include:</u> Material puncture/hole/rupture, sharp edges, material twisted/bent, difficult to remove, misassembled during manufacturing, <u>Common injuries resulted from:</u> Burst container/vessel, difficult to remove, material puncture/hole/rupture, device handling, material fragmentation



Review of Post Market Surveillance / Customer Feedback Data

Post-Market Surveillance Data	Summary of activity/results
Customer complaints	4 product complaints reported. 2 complaints confirmed to involve Silq catheter defects (details in Review of Product Complaints section) • Agglomeration on catheters • Retaining sleeve fell off catheter
Customer feedback	None
Product returns	None
Field Safety Corrective Action/Recalls	None
Public health advisories of safety notices from regulatory authorities	None
Nonconforming product	None
Device-related CAPA's	None



Review of Post Market Surveillance / Customer Feedback Data

Post-Market Surveillance Data	Summary of activity/results
Incoming, in-process and final quality inspections	0 NCs found at receiving inspection 0 NCs found at in-process inspection 0 NCs found at final inspection 0 NCs found in field 0 NCs found in inventory / finished goods
Manufacturing yields	<u>Suspension</u> None <u>Finished device:</u> 3 of 4 production lots accepted; 60.3 cases rejected (72% yield) <u>2-Way 14Fr x 10mL, Sterile, 10pack: Released 22.1 cases</u> <u>2-Way 16Fr x 10mL, Sterile, 10pack: Released 94.4 cases</u> <u>2-Way 18Fr x 10mL, Sterile, 10pack: Released 38 cases</u>



Review of Post Market Surveillance / Customer Feedback Data

Post-Market Surveillance Data	Summary of activity/results
Process validations	Coating 14Fr and 18Fr process validation on automated coating machine completed in 2023. Report released in May 2024 (TR-SLQ-014) Sealing process with Impulse Sealer PQ completed with another pouch manufacturer in 2023. Report released in 2024 (TR-SLQ-0011)
Supplier issues	None
Training issues	None



Review of reportings to regulatory authorities



Review of Reportings to Regulatory Authorities

- There have been no reportings to regulatory authorities in 2024.



Review of Statistical Techniques



Review of Statistical Techniques

- Statistical techniques is documented in QM.SLQ011 Rev A.
 - Suitable and adequate for its application to:
 - Design Verification & Validation
 - Process qualifications and validation
- Suspension Production sampling plans (Silq)
 - 100% QC analysis per QC-C.SLQ001 (appearance, residue on evaporation, pH, ^1H NMR, ^{19}F NMR)



Review of Statistical Techniques

- Finished device inspection sampling plans (Pathway)
 - Sampling based on C=0 Sampling Plan. No single defective unit per designated sample size. (QP-030 Rev D Statistical Techniques)
 - AQL assigned based on severity rating

Table 7.5.3 – AQL Determination		
AQL	Risk Analysis Severity Rankings	
	Severity Rating	Definition
6.5	1	Negligible – No injury to patient or user.
4.0	2	Minor – Temporary injury or impairment not requiring professional medical intervention.
2.5	3	Moderate – Moderate injury that may require medical intervention.
1.0	4	Serious – Serious injury requiring medical attention to prevent permanent impairment.
0.65/ 0.25	5	Critical – Life threatening injury, disability, death



Review of supplier performance



Review of Supplier Performance

Category I suppliers

- Pathway (Approved): Re-evaluation in-process
- Richman Chemical (Approved): NO PO's were issued to Richman Chemical in 2024. Supplier re-evaluation was performed using previous re-evaluation score from 2023. No change to their approval status.
- TCI (Approved): 1 PO was issued to TCI in 2024. The PO is currently in-process. Supplier re-evaluation was performed using previous re-evaluation score from 2023. No change to their approval status.
- HPFY (Approved): NO PO's were issued to HPFY in 2024. Supplier re-evaluation was performed using the initial self-assessment survey. No change to their approval status.
- Steripax (Approved): NO PO's were issued to Steripax in 2024. Supplier re-evaluation was performed using the initial self-assessment survey. No change to their approval status.



Review of Supplier Performance

Category III suppliers

- Fisher Scientific (Approved): No PO's were issued to Fisher Scientific during the evaluation period (Nov 2022 - Dec 2024). Supplier re-evaluation was performed using the initial supplier assessment score. No change to their approval status.
- Independent Air Group (Approved): No PO's were issued to Independent Air Group during the evaluation period (Dec 2022 - Dec 2024). Supplier re-evaluation was performed using the initial supplier assessment score. No change to their approval status.
- McMaster-Carr (Approved): No PO's were issued to McMaster-Carr during the evaluation period (Nov 2022 - Dec 2024). Supplier re-evaluation was performed using the initial supplier assessment score. No change to their approval status.
- Micro Precision Calibration (Approved): 5 PO's were issued to MPC during the evaluation period (Aug 2022 - Dec 2024). Supplier re-evaluation was performed. Scored 100/100. No change to their approval status.
- Repligen (Approved): 1 PO was issued to Repligen during the evaluation period (Aug 2022 – Dec 2024). Supplier re-evaluation was performed. Scored 99.5/100. No change to their approval status.



Review of Supplier Performance

Category III suppliers (cont.)

- VWR (Approved): No PO's were issued to VWR during the evaluation period (Nov 2022 - Dec 2024). Supplier re-evaluation was performed using the initial supplier assessment score. No change to their approval status.



Review of Supplier Performance

Category IV suppliers

- Uline (Approved): No PO's were issued to Uline during the evaluation period (Nov 2022 - Dec 2024). Supplier re-evaluation was performed using the initial supplier assessment score. No change to their approval status.
- United States Plastic Corp (Approved): 1 PO was issued to US Plastic Corp during the evaluation period (Nov 2022 – Dec 2024). Supplier re-evaluation was performed. Scored 100/100. No change to their approval status.
- Cole Parmer (Approved): 1 PO was issued to Cole Parmer during the evaluation period (Nov 2022 – Dec 2024). Supplier re-evaluation was performed. Scored 100/100. No change to their approval status.
- Glacier Tanks (Approved): 1 PO was issued to Glacier Tanks during the evaluation period (Nov 2022 – Dec 2024). Supplier re-evaluation was performed. Scored 100/100. No change to their approval status.



Changes that could affect the QMS



Changes that could affect QMS

- New or revised regulatory standards
 - ISO 15223-1:2021 replaces ISO 15223-1:2016 (Medical devices — Symbols to be used with information to be supplied by the manufacturer, Part 1: General requirements)
 - Addition of symbol for UDI (Unique Device Identifier)
 - Product labels for pouch and box of Silq product (SLQ-LB-4000, SLQ-LB-4002, SLQ-LB-4008, SLQ-LB-4009) contain UDI. No action required.
 - ASTM D4169-23e1 replaces ASTM D4169-22 (Standard Practice for Performance Testing of Shipping Containers and Systems)
 - Changes to Schedule C (Vehicle Stacking) for small packages (smaller than 2.0 ft³ and weighs between 10.0 lbs and 30.0 lbs). Does not apply to the corrugated container used in the packaging validation testing (1.25 ft³, 7.88lb).
 - No impact on testing performed in accordance with old revision of standard
 - ASTM D4332-22 replaces ASTM D41332-14 (Standard Practice for Conditioning Containers, Packages, or Packaging Components for Testing)
 - No significant changes
 - No impact on testing performed in accordance with old revision of standard
 - ASTM F1886/F1886M-16(2024) replaces ASTM F1886/F1886M-16 (Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection)
 - No changes. 2024 revision is a reapproval of the 2016 revision.
 - No impact on testing performed in accordance with old revision of standard



Changes that could affect QMS

- Recommendations for improvement
 - Product
 - Complete Gage R&R study on the following in-process quality control test methods: UV Light Absorbance (TM-HDX-001) and Rinse Test (TM-HDX-002).
 - Complete test method validation for ^1H and ^{19}F NMR spectra readings



Review of Risk Management Activities



Review of Risk Management Activities

- Based on review of data obtained from post-market surveillance activities, no new hazards have been identified. Potential change to an existing hazard #10.1 regarding three failure modes:
 1. Balloon deflates or does not inflate – Valve malfunction or leaks (wall thickness issue)
 2. Balloon does not inflate – valve malfunctions
 3. Catheter breaks under normal tensile stresses – Material or dimensional issues
- The probability of each failure mode is currently set at '1'. *Action item to evaluate the probability rating and determine if the rating should increase based on complaints seen in the field.*



Review of Design Control Activities



Review of Design Control Activities

Shelf-life extension for Suspension (C.SLQ001)

- Study performed in accordance with Protocol #VV.SLQ028
- Catheter samples coated with 2 year real-time aged Suspension C.SLQ001 (Lot #HXY-825021)
- Coating adherence verification and crystal violet absorbance test were performed on coated samples. Results met acceptance criteria confirming a shelf life of 2 years for Suspension. Report #VV.029.



Review of Resource Requirements



Review of Resource Requirements

- Silq
 - No change in the organization/personnel in 2024.
 - Transition quality and regulatory responsibilities to Pathway in 2025.
- CMO's
 - Pathway: Manufacture of coated catheters, labeling, final packaging, sterilization and distribution.



End of Management Review

Next Management Review will be held in January 2026

